

THE TRUSTED NAME IN HEALTH SUPPORTS FOR 125 YEARS

SAI  **THERAPEUTIC
BRANDS**

SUPPORTS AND BRACES



**OTC PROFESSIONAL ORTHOPAEDIC
CHAMPION HEALTH AND SPORTS SUPPORTS**

SAI THERAPEUTIC BRANDS

THE SAI THERAPEUTIC BRANDS TRADITION



For well over a century, SAI Therapeutic Brands has provided premium quality, cost effective healthcare products to professional orthotic fitters worldwide. Throughout, we have labored intensely to build an industry-wide reputation where quality in materials and workmanship, and dealer education are synonymous with our TRUFORM® Compression Hosiery, OTC® Professional Orthopaedic, CHAMPION® Health and Sports Supports, PCP® Home Healthcare & Mobility, CSX Athletic Braces, Supports, & Compression Garments, and TOUCH Fashionable Compression Socks brands.

While our history is rich with new product development and innovation, reinvention has long been the real key to our continued success, and has enabled us to maintain our long-time role as a full-line industry leader in the design, manufacture and distribution of orthopaedic soft goods, prescription grade bracing, and rigid motion control orthoses.

Today the tradition of reinvention and innovation continues as we continue to refine and improve our complement of handcrafted appliances: The OTC® Professional Orthopaedic and CHAMPION® lines of consumer-friendly supports and braces!



VALUE

You don't have to pay more to get the very best! High quality and affordability are synonymous with OTC and Champion products - a real plus in today's competitive market.



CONSISTENCY

In design and quality - with few exceptions you buy direct from the original manufacturer, which ensures consistent product quality from order to order.



SIGNATURE CUSTOMER SERVICE

95% of all orders are shipped within 24 hours of receipt, with a 99% accuracy rate - you receive what you want, when you want!



SPECIAL ORDER CAPABILITY

Orthopedic appliances can be special-ordered to fit unique figure types.



WARRANTY

OTC and Champion products are warranted to be free of material or workmanship defects. Problematic merchandise will be repaired or replaced.



FULL LINE SUPPLIER - 1 STOP SHOP

Virtually every imaginable client need can be satisfied from the SAI Brands program! Literally the best mix of motion control, support, safety convenience and mobility products, all from one source of supply!



DEDICATED SALES ASSOCIATES

Fully trained professionals will personally visit your facility on a regular basis to assist in managing various aspects of your business.

CONTENTS:



Knee Supports	4
Wrist Supports	5
Cervical and Clavicle Supports	5
Rib Supports	6
Abdominal and Sacroiliac Supports	6



Knee Supports	7
Ankle and Foot Supports	9
Elbow Straps and Supports	12
Wrist Supports	13
Maternity and Back Supports	15
Rib Supports	17
Abdominal Hernia Supports	17
Cervical Clavicle Supports	18
Arthritic Warmers	21



Knee Supports	21
Ankle Supports	25
Elbow Straps and Supports	26
Wrist Supports	27
Shoulder Supports	29
Back and Rib Supports	30
Hernia and Sport Supports	31
Hot/Cold Compresses	31
Bandages and Cast Covers	32
Pediatric Supports	33



SAI Tradition	1
Planograms	34
Bracing and Supports Product Index	35



Throughout this catalog, the level of support that each product provides is shown to assist in choosing the right item for individual patient need as follows:

- **MINIMUM** - Provides compressive support without any reinforcement. For weak and aching muscles and joints, use after recovery or to prevent injury.
- **MEDIUM** - Provides compressive support along with semi-rigid reinforcement. For mild to moderate conditions that might require some stabilization. For post-injury use.
- **MAXIMUM** - Provides compressive support along with rigid reinforcement. For moderate to severe conditions requiring motion control to the affected area. For immediate post-injury use.

PRACTITIONER AND PATIENT FOCUSED

The SAI Therapeutic Brands program is built on the best mix of products that are designed and engineered for therapeutic performance and patient comfort - and specifically targeted to meet the needs of the everyday active consumer; the younger, fitness-oriented consumer with self care/cash-pay orientation; and the occupational injury consumer.

Dynamic contemporary packaging educates the consumer, and creates a bold retail presence to help generate faster turns. Innovative packaging design features create planograms that are easy to shop and easy to merchandise.



Size clearly labeled

Product photo clearly shows all features and components

Consumer bullet copy ties product features together with user end benefits

Care and application instructions on side panel

Degree of support



Product features

Spanish language on side panel with features and benefits information, care and application instructions

Typical indications for product

QR Code to access product demonstration video

Where to measure product and how to choose correct size



KNEE SUPPORTS

OTC SELECT IS THE RESULT OF EXTENSIVE RESEARCH into the latest high-tech materials and fabrication processes. All products feature an ideal combination of wicking, moisture management and odor reducing anti-microbial treatments to ensure all-day wearing comfort.

OTC SELECT PRODUCTS INCLUDE:

- **ANTI-MICROBIAL, ANTI-BACTERIAL TREATED LATEX-FREE MATERIALS** that maintain breathability and freshness, and fight odor even after prolonged periods of wear.
- **NON-ABRASIVE HOOK-AND-LOOP FASTENING** materials with no sharp edges to irritate skin or damage clothing, yet provide for superior gripping over thousands of applications.
- **HIGHLY CONTROLLED SEWING PROCESSES** that result in tighter stitching, lower profile seams and a sleek, contemporary look.
- **BLACK AND SLATE COLORS** for a less clinical, more contemporary appearance.

MODEL: 2541



ORTHOTEX KNEE STABILIZER WITH SPIRAL STAYS

HPCPS: A4466 •• MEDIUM

- Spiral stays provide moderate medial/lateral stability
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2541-S	S	12.5-13.75 in / 31.7-34.9 cm
2541-M	M	14-15.25 in / 35.5-38.7 cm
2541-L	L	15.5-16.5 in / 39.3-41.9 cm
2541-XL	XL	16-17.5 in / 40.6-44.4 cm
2541-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
2541-3XL	3XL	19.5-21 in / 49.5-53.3 cm
2541-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 2542



ORTHOTEX KNEE STABILIZER WRAP WITH SPIRAL STAYS

HPCPS: A4466 •• MEDIUM

- Opens flat for easy application over injured, swollen or tender knees
- Spiral stays provide moderate medial/lateral stability
- Encircling straps are reversible for wear on right or left knee

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2542-XS	XS	12-13 in / 20.5-33 cm
2542-S	S	13-14 in / 33-35.5 cm
2542-M	M	14-15 in / 35.5-38.1 cm
2542-L	L	15-16 in / 38.1-40.6 cm
2542-XL	XL	16-17 in / 40.6-43.2 cm
2542-2XL	2XL	17-18 in / 43.2-45.7 cm
2542-3XL	3XL	18-19 in / 45.7-48.3 cm
2542-4XL	4XL	19-20 in / 48.3-50.8 cm
2542-5XL	5XL	20-22 in / 50.8-55.9 cm

MODEL: 2543



ORTHOTEX KNEE STABILIZER WITH HINGED BARS

HPCPS: L1812 ••• MAXIMUM

- Hinged bars provide optimum medial/lateral stability
- Multifunction buttress surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee

INDICATIONS: Medial/lateral knee instability, Chondromalacia, osteoarthritis, compression and support

MEASURE AROUND THE BEND OF THE KNEE

2543-S	S	12.5-13.75 in / 31.7-34.9 cm
2543-M	M	14-15.25 in / 35.5-38.7 cm
2543-L	L	15.5-16.5 in / 39.3-41.9 cm
2543-XL	XL	16-17.5 in / 40.6-44.4 cm
2543-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
2543-3XL	3XL	19.5-21 in / 49.5-53.3 cm
2543-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 2544



ORTHOTEX KNEE STABILIZER WRAP WITH HINGED BARS

HPCPS: L1812 ••• MAXIMUM

- Opens flat for easy application over injured, swollen or tender knees
- Hinged bars provide optimum medial/lateral stability
- Encircling straps are reversible for wear on right or left knee

INDICATIONS: Medial/lateral knee instability, Chondromalacia, osteoarthritis, compression and support

MEASURE AROUND THE BEND OF THE KNEE

2544-XS	XS	12-13 in / 20.5-33 cm
2544-S	S	13-14 in / 33-35.5 cm
2544-M	M	14-15 in / 35.5-38.1 cm
2544-L	L	15-16 in / 38.1-40.6 cm
2544-XL	XL	16-17 in / 40.6-43.2 cm
2544-2XL	2XL	17-18 in / 43.2-45.7 cm
2544-3XL	3XL	18-19 in / 45.7-48.3 cm
2544-4XL	4XL	19-20 in / 48.3-50.8 cm
2544-5XL	5XL	20-22 in / 50.8-55.9 cm

MODEL: 2546



ORTHOTEX KNEE SUPPORT WITH STABILIZER PAD

HPCPS: A4466 • MINIMUM

- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Multi-directional stretch materials provide even compression over surface anatomy

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2546-S	S	12.5-13.75 in / 31.7-34.9 cm
2546-M	M	14-15.25 in / 35.5-38.7 cm
2546-L	L	15.5-16.5 in / 39.3-41.9 cm
2546-XL	XL	16-17.5 in / 40.6-44.4 cm
2546-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
2546-3XL	3XL	19.5-21 in / 49.5-53.3 cm
2546-4XL	4XL	21-24 in / 53.3-61 cm



WRIST AND THUMB SUPPORTS



MODEL: 2382



Specify for RIGHT or LEFT wrist

6' WRIST SPLINT

HCPCS: L3908 ••• MAXIMUM

- For petite to average females and adolescents
- Lightweight glove with separate cinch fasteners provides customized support and compression
- Low contour at palmar crease allows for full finger dexterity and grip strength
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2382-XS	XS	4.5-5.5 in / 11.4-14 cm
2382-S	S	5.5-6.5 in / 14-16.5 cm
2382-M	M	6.5-7.5 in / 16.5-19 cm
2382-L	L	7.5-8.5 in / 19-21.6 cm
2382-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2386



Specify for RIGHT or LEFT wrist

6' WRIST-THUMB SPLINT

HCPCS: L3908 ••• MAXIMUM

- For petite to average females and adolescents
- Lightweight glove with separate cinch fasteners provides customized support and compression
- Radial splint anchored by proximal circumference strap provides for optimum immobilization of the thumb
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: De Quervain's tendonitis, strains and sprains of the MCP joint, basal joint arthritis, advanced carpal tunnel syndrome, cumulative trauma disorders of the thumb and wrist

MEASURE AROUND THE BREAK OF THE WRIST

2386-XS	XS	4.5-5.5 in / 11.4-14 cm
2386-S	S	5.5-6.5 in / 14-16.5 cm
2386-M	M	6.5-7.5 in / 16.5-19 cm
2386-L	L	7.5-8.5 in / 19-21.6 cm
2386-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2383



Specify for RIGHT or LEFT wrist

8' WRIST SPLINT

HCPCS: L3908 ••• MAXIMUM

- For average to larger adults
- Lightweight glove with separate cinch fasteners provides customized support and compression
- Low contour at palmar crease allows for full finger dexterity and grip strength
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2383-XS	XS	4.5-5.5 in / 11.4-14 cm
2383-S	S	5.5-6.5 in / 14-16.5 cm
2383-M	M	6.5-7.5 in / 16.5-19 cm
2383-L	L	7.5-8.5 in / 19-21.6 cm
2383-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2387



Specify for RIGHT or LEFT wrist

8' WRIST-THUMB SPLINT

HCPCS: L3908 ••• MAXIMUM

- For average to larger adults
- Lightweight glove with separate cinch fasteners provides customized support and compression
- Radial splint anchored by proximal circumference strap provides for optimum immobilization of the thumb
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: De Quervain's tendonitis, strains and sprains of the MCP joint, basal joint arthritis, advanced carpal tunnel syndrome, cumulative trauma disorders of the thumb and wrist

MEASURE AROUND THE BREAK OF THE WRIST

2387-XS	XS	4.5-5.5 in / 11.4-14 cm
2387-S	S	5.5-6.5 in / 14-16.5 cm
2387-M	M	6.5-7.5 in / 16.5-19 cm
2387-L	L	7.5-8.5 in / 19-21.6 cm
2387-XL	XL	8.5-10 in / 21.6-25.4 cm

CERVICAL AND CLAVICLE SUPPORTS



MODEL: 2454



FIGURE 8 CLAVICLE STRAP

HCPCS: L3650 •• MEDIUM

- Holds the shoulder region in the correct position to help promote healing following a clavicle fracture or separation
- Figure 8 design for easy application
- Varying degrees of support can be obtained by adjusting the tension of the shoulder straps

INDICATIONS: Clavicle fracture or separation, poor posture

MEASURE AROUND THE CHEST AT THE LEVEL OF THE UNDERARMS

2454-XS	XS	16-24 in / 50.8-61 cm
2454-S	S	24-30 in / 61-76.2 cm
2454-M	M	30-36 in / 76.2-91.4 cm
2454-L	L	36-42 in / 91.4-106.7 cm
2454-XL	XL	42-48 in / 106.7-121.9 cm

MODEL: 2503



INFLATABLE CERVICAL TRACTION UNIT

HCPCS: E0856 ••• MAXIMUM

- Easy to use, comfortable and portable
- Unit wraps around the neck; air chambers are inflated with easy-to-use hand pump, elongating to relieve pressure over cervical muscles and nerve tissue

INDICATIONS: Musculoskeletal or neurological impairment in the cervical region, osteoarthritis, cervical spondylosis

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2503	UNIVERSAL
------	-----------





RIB SUPPORTS

MODEL: 2458



RIB BELT FOR MEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Ribbed web elastic molds around surface anatomy for a comfortable, custom fit
- 6" deep all around

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2458-UR	UR	29 – 38 in / 73.6-96.5 cm
2458-UL	UL	39 – 48 in / 99-121.9 cm

MODEL: 2658



RIB BELT FOR WOMEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Ribbed web elastic molds around surface anatomy for a comfortable, custom fit
- Cut low in front to accommodate the bust line
- 4" in front, 6" in back

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2658-UR	UR	29 – 38 in / 73.6-96.5 cm
2658-UL	UL	39 – 48 in / 99-121.9 cm



ABDOMINAL AND SACROILIAC SUPPORTS

MODEL: 2508



MULTIPLE USE ELASTIC BINDER FOR MEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Ribbed web elastic molds around surface anatomy for a comfortable, custom fit
- Anti-microbial fabric, non-abrasive fastening
- Can also be used as a rib fracture support
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2508-XS	XS	25-28 in / 64-71 cm
2508-S	S	29-32 in / 73.6-81.2 cm
2508-M	M	33-36 in / 83.8-91.4 cm
2508-L	L	37-40 in / 93.9-101.6 cm
2508-XL	XL	41-44 in / 104.1-111.8 cm
2508-2XL	2XL	45-48 in / 114.3-121.9 cm

MODEL: 2688



MULTIPLE USE ELASTIC BINDER FOR WOMEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Ribbed web elastic molds around surface anatomy for a comfortable, custom fit
- Anti-microbial fabric, non-abrasive fastening
- Contoured to fit the female form
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2688-S	S	19-22 in / 48.3-55.9 cm
2688-M	M	23-26 in / 58.4-66 cm
2688-L	L	27-30 in / 68.6-76.2 cm
2688-XL	XL	31-34 in / 78.7-86.4 cm
2688-2XL	2XL	35-38 in / 88.9-96.5 cm

MODEL: 2955



ABDOMINAL HERNIA BELT

HPCPS: A4466 • MINIMUM

- Ribbed web elastic molds around surface anatomy for a comfortable, custom fit
- Comfortable molded pad provides localized support at the umbilicus
- Anti-microbial fabric, non-abrasive fastening
- 6" deep all around

INDICATIONS: Umbilical hernia, weak or pendulous abdomen, abdominal contusions and tenderness

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2955-XS	XS	25-28 in / 64-71 cm
2955-S	S	29-32 in / 73.6-81.2 cm
2955-M	M	33-36 in / 83.8-91.4 cm
2955-L	L	37-40 in / 93.9-101.6 cm
2955-XL	XL	41-44 in / 104.1-111.8 cm
2955-2XL	2XL	45-48 in / 114.3-121.9 cm

MODEL: 2885



SACROILIAC SUPPORT / TROCHANTER BELT

HPCPS: A4466 • MINIMUM

- Provides firm support and localized circular compression
- Bilateral tightening pull straps with non-abrasive fastening
- Easy to apply, comfortable to wear

INDICATIONS: Support and compression for the sacroiliac joint, arthritis, joint hypermobility, inequality of leg length, maternity support

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2885	UNIVERSAL
------	-----------



KNEE SUPPORTS



MODEL: 2416



PULLOVER ELASTIC KNEE SUPPORT

HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Anatomic design works in concert with the motion of the knee, encircles like a bandage but without the bulk

INDICATIONS: Mild sprains and strains, joint weakness due to previous injury or overexertion

MEASURE AROUND THE BEND OF THE KNEE

2416-S	S	10-12.75 in / 25.4-32.4 cm
2416-M	M	13-15.75 in / 33-40 cm
2416-L	L	16-18.75 in / 40.6-47.6 cm
2416-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 2422



KNEEDIT® THERAPEUTIC KNEE GUARD

HCPCS: A4466 •• MEDIUM

- Patented anatomic design provides support, comfort and mobility during activity
- Semi-rigid support with high density padding molds around soft tissues for customized support
- Superior therapy that can be used in conjunction with traditional methods like taping and wrapping

INDICATIONS: Arthritis, tendonitis, Chondromalacia, knee strains

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2422	UNIVERSAL
------	-----------

MODEL: 2425



KNEE SUPPORT WITH VISCOELASTIC INSERT

HCPCS: A4466 •• MEDIUM

- Anatomically sculpted insert pad provides precise compression to support weak or injured knees
- Controlled, four-way stretch elastic offers effective compressive support to the entire knee area
- Spiral stays provide moderate medial/lateral stability

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2425-S	S	14-16 in / 35.6-40.6 cm
2425-M	M	16-18 in / 40.6-45.7 cm
2425-L	L	18-20 in / 45.7-50.8 cm
2425-XL	XL	20-22 in / 50.8-55.9 cm

MODEL: 2548



ORTHOTEX KNEE STABILIZER WITH RANGE OF MOTION HINGED BARS

HCPCS: L1833 ••• MAXIMUM

- Low profile single axis hinged bars provide medial-lateral stability, adjust to control knee extension at 15°, 30° or 50°
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization

INDICATIONS: Post-operative rehabilitation, patellar malalignment and instability, cruciate and collateral ligament support, medial/lateral cartilage support, patellar tendonitis

MEASURE AROUND THE BEND OF THE KNEE

2548-S	S	12.5-13.75 in / 31.7-34.9 cm
2548-M	M	14-15.25 in / 35.5-38.7 cm
2548-L	L	15.5-16.5 in / 39.3-41.9 cm
2548-XL	XL	16-17.5 in / 40.6-44.4 cm
2548-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
2548-3XL	3XL	19.5-21 in / 49.5-53.3 cm
2548-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 2540



ORTHOTEX KNEE STABILIZER WRAP FOR OA WITH HINGED BARS

HCPCS: 1820 ••• MAXIMUM

- Elastic wraparound strap system with single polycentric hinged bar provides counter resistance and exerts corrective pressure to offload affected side
- Interior condyle pads with different thicknesses can be adjusted for varying daily needs

INDICATIONS: Osteoarthritis, lateral pressure syndrome, patella instabilities, lateral patella malalignment

MEASURE AROUND THE BEND OF THE KNEE

2540-XS	XS	12-13 in / 20.5-33 cm
2540-S	S	13-14 in / 33-35.5 cm
2540-M	M	14-15 in / 35.5-38.1 cm
2540-L	L	15-16 in / 38.1-40.6 cm
2540-XL	XL	16-17 in / 40.6-43.2 cm
2540-2XL	2XL	17-18 in / 43.2-45.7 cm
2540-3XL	3XL	18-19 in / 45.7-48.3 cm
2540-4XL	4XL	19-20 in / 48.3-50.8 cm
2540-5XL	5XL	20-22 in / 50.8-55.9 cm

Specify for RIGHT or LEFT knee

MODEL: 2549



ORTHOTEX STABILIZER WRAP WITH RANGE OF MOTION HINGED BARS

HCPCS: 1833 ••• MAXIMUM

- Opens flat for easy application over injured, swollen or tender knees
- Low profile single axis hinged bars provide medial-lateral stability, adjust to control knee extension at 15°, 30° or 50°
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization

INDICATIONS: Post-operative rehabilitation, patellar malalignment and instability, cruciate and collateral ligament support, medial/lateral cartilage support, patellar tendonitis

MEASURE AROUND THE BEND OF THE KNEE

2549-XS	XS	12-13 in / 20.5-33 cm
2549-S	S	13-14 in / 33-35.5 cm
2549-M	M	14-15 in / 35.5-38.1 cm
2549-L	L	15-16 in / 38.1-40.6 cm
2549-XL	XL	16-17 in / 40.6-43.2 cm
2549-2XL	2XL	17-18 in / 43.2-45.7 cm
2549-3XL	3XL	18-19 in / 45.7-48.3 cm
2549-4XL	4XL	19-20 in / 48.3-50.8 cm
2549-5XL	5XL	20-22 in / 50.8-55.9 cm





KNEE SUPPORTS

MODEL: 2545



KNEE SUPPORT WITH FRONT OPENING HCPCS: A4466 •• MEDIUM

- Opens flat for easy application over swollen or tender knees
- Strong controlled stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Bilateral pressure pads cradle the patella for added stability and to help maintain proper patella alignment

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2545-XS	XS	9.5-11.5 in / 24.1-29.2 cm
2545-S	S	11.5-13 in / 29.2-33 cm
2545-M	M	13-14.75 in / 33-37.5 cm
2545-L	L	14.75-16.5 in / 37.5-41.9 cm
2545-XL	XL	16.5-18.5 in / 41.9-47 cm
2545-2XL	2XL	18.5-21 in / 47-53.3 cm
2545-3XL	3XL	21-23.5 in / 53.3-59.7 cm

MODEL: 2553



KNEE SUPPORT WITH SPIRAL STAYS HCPCS: A4466 •• MEDIUM

- Spiral stays provide moderate medial/lateral stability
- Strong controlled stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Excellent for post-surgical recovery, increases wearer confidence for walking or standing

INDICATIONS: Sprains and strains, osteoarthritis, general support

MEASURE AROUND THE BEND OF THE KNEE

2553-XS	XS	9.5-11.5 in / 24.1-29.2 cm
2553-S	S	11.5-13 in / 29.2-33 cm
2553-M	M	13-14.75 in / 33-37.5 cm
2553-L	L	14.75-16.5 in / 37.5-41.9 cm
2553-XL	XL	16.5-18.5 in / 41.9-47 cm
2553-2XL	2XL	18.5-21 in / 47-53.3 cm
2553-3XL	3XL	21-23.5 in / 53.3-59.7 cm

MODEL: 2554



KNEE SUPPORT WITH HINGED BARS HCPCS: L1812 ••• MAXIMUM

- Hinged bars for optimum medial/lateral stability
- Strong controlled stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Provides protection and stability for post-surgical recovery, or for those with a history of joint instability

INDICATIONS: Medial/lateral knee instability, compression and support

MEASURE AROUND THE BEND OF THE KNEE

2554-XS	XS	9.5-11.5 in / 24.1-29.2 cm
2554-S	S	11.5-13 in / 29.2-33 cm
2554-M	M	13-14.75 in / 33-37.5 cm
2554-L	L	14.75-16.5 in / 37.5-41.9 cm
2554-XL	XL	16.5-18.5 in / 41.9-47 cm
2554-2XL	2XL	18.5-21 in / 47-53.3 cm
2554-3XL	3XL	21-23.5 in / 53.3-59.7 cm

MODEL: 2555



KNEE SUPPORT WITH PRESSURE PADS HCPCS: A4466 •• MEDIUM

- Bilateral pressure pads provide compression and stabilization around patella
- Strong controlled stretch elastic for effective compressive support over the knee along with all-day wearing comfort
- An excellent support for post-surgical recovery, it helps reduce the incidence of reinjury while active

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2555-XS	XS	9.5-11.5 in / 24.1-29.2 cm
2555-S	S	11.5-13 in / 29.2-33 cm
2555-M	M	13-14.75 in / 33-37.5 cm
2555-L	L	14.75-16.5 in / 37.5-41.9 cm
2555-XL	XL	16.5-18.5 in / 41.9-47 cm
2555-2XL	2XL	18.5-21 in / 47-53.3 cm
2555-3XL	3XL	21-23.5 in / 53.3-59.7 cm

MODEL: 2555V



KNEE SUPPORT WITH PRESSURE PADS AND EXPANSION PANEL HCPCS: A4466 •• MEDIUM

- Gusseted opening with opposing hook-and-loop fastening straps permit easy application for those with a well developed, full size thigh
- Strong controlled stretch elastic for effective compressive support over the knee along with all-day wearing comfort
- Bilateral pressure pads provide compression and stabilization around patella

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

2555V-XS	XS	9.5-11.5 in / 24.1-29.2 cm
2555V-S	S	11.5-13 in / 29.2-33 cm
2555V-M	M	13-14.75 in / 33-37.5 cm
2555V-L	L	14.75-16.5 in / 37.5-41.9 cm
2555V-XL	XL	16.5-18.5 in / 41.9-47 cm
2555V-2XL	2XL	18.5-21 in / 47-53.3 cm
2555V-3XL	3XL	21-23.5 in / 53.3-59.7 cm



ANKLE AND FOOT SUPPORTS



MODEL: 2092



ANKLE STIRRUP BRACE WITH AIRFORM PADS

HPCPS: L4350 ••• MAXIMUM

- Unique multi-chamber air bladder captures and provides a smooth interface with the surface anatomy around the ankle
- Rigid plastic shells provide protection and help prevent inversion and eversion
- Quick adjust heel loops allow for a custom fit and permit unobstructed plantar- and dorsi-flexion.

INDICATIONS: Acute ankle injuries, tender or swollen ankles

UNIVERSAL FIT

2092	UNIVERSAL
2092-P	YOUTH (PONY)

Specify for RIGHT or LEFT ankle

MODEL: 2371



HIGH PERFORMANCE ANKLE BRACE

HPCPS: L1902 •• MEDIUM

- Two bilateral spiral stays provide moderate medial/lateral stability
- Strong, lightweight flannel lined canvas for a firm, comfortable fit
- Elastic panel in back allows for normal athletic foot movement like running and jumping

INDICATIONS: Acute ankle sprains, chronic ankle instability

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2371-XS	XS	5.25-6.5 in / 13.3-16.5 cm
2371-S	S	6.5-7.75 in / 16.5-19.7 cm
2371-M	M	7.75-9 in / 19.7-22.9 cm
2371-L	L	9-10.25 in / 22.9-26 cm
2371-XL	XL	10.25-11 in / 26-27.9 cm

MODEL: 2372



CANVAS ANKLE SPLINT

HPCPS: L1902 ••• MAXIMUM

- Rigid pre-bent metal stays are shapeable and provide superior medial-lateral stability
- Strong, lightweight flannel lined canvas for a firm, comfortable fit
- Effective in stabilizing the ankle following severe sprains, tendonitis or after cast removal

INDICATIONS: Acute ankle sprains, chronic ankle instability, post cast removal

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2372-XS	XS	5.25-6.5 in / 13.3-16.5 cm
2372-S	S	6.5-7.75 in / 16.5-19.7 cm
2372-M	M	7.75-9 in / 19.7-22.9 cm
2372-L	L	9-10.25 in / 22.9-26 cm
2372-XL	XL	10.25-11 in / 26-27.9 cm

MODEL: 2375



ANKLE STABILIZER WITH HEEL LOCKING STRAP

HPCPS: L1902 ••• MAXIMUM

- Adjustable strapping system provides optimum medial-lateral control, and will not stretch during exercise
- Lightweight materials allow full range of dorsi- and plantar-flexion
- Can be worn after an acute ankle sprain or as a protection against repeat ankle sprains during athletic activities

INDICATIONS: Acute ankle sprains, chronic ankle instability, post cast removal

MEASURE AROUND THE ANKLE, FROM THE BACK/BOTTOM OF THE HEEL AROUND THE FRONT OF THE ANKLE

2375-XS	XS	10-11 in / 25.4-27.9 cm
2375-S	S	11-12 in / 27.9-30.5 cm
2375-M	M	12-13 in / 30.5-33 cm
2375-L	L	13-14 in / 33-35.5 cm
2375-XL	XL	14-15 in / 35.5-38.1 cm
2375-2XL	2XL	15-17 in / 38.1-43.2 cm

MODEL: 2376



ANKLE STABILIZER WITH EXOSKELETON AND STRAPS

HPCPS: L1902 ••• MAXIMUM

- Combines the compression and comfort of a soft ankle brace with the rigid stability of an ankle stirrup
- Plastic exoskeleton is formed directly onto the soft breathable fabric, forming a network of strategically placed support zones
- Fully adjustable strapping system provides optimum medial-lateral control, will not stretch during exercise

INDICATIONS: Acute ankle sprains, chronic ankle instability, post cast removal

SHOE SIZE

		MEN'S	WOMEN'S
2376-XS	XS	5.5 - 7	6.5 - 8
2376-S	S	7.5 - 9	8.5 - 10
2376-M	M	9.5 - 11	10.5 - 12
2376-L	L	11.5 - 13	12.5 - 14
2376-XL	XL	13.5 - 14	N/A

MODEL: 2417



PULLOVER ELASTIC ANKLE SUPPORT

HPCPS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the ankle along with all-day wearing comfort
- Anatomic design works in concert with the motion of the ankle, encircles like a bandage but without the bulk

INDICATIONS: Aching joints and muscles, edema, joint tenderness

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2417-S	S	6.75-7.75 in / 17.1-19.7 cm
2417-M	M	7.75-8.75 in / 19.7-22.2 cm
2417-L	L	8.75-9.75 in / 22.2-24.8 cm
2417-XL	XL	9.75-10.75 in / 24.8-27.3 cm





ANKLE AND FOOT SUPPORTS

MODEL: 2426



ANKLE SUPPORT WITH VISCOELASTIC INSERT

HPCPS: A4466 ●● MEDIUM

- Anatomically sculpted viscoelastic silicone insert delivers target support to weakened or injured soft tissues around the ankle
- Controlled, four-way stretch elastic offers effective compressive support to the entire knee area

INDICATIONS: Swollen or tender ankles resulting from sprains, strains or sports injuries

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2426-S	S	6.75-7.75 in / 17.1-19.7 cm
2426-M	M	7.75-8.75 in / 19.7-22.2 cm
2426-L	L	8.75-9.75 in / 22.2-24.8 cm
2426-XL	XL	9.75-10.75 in / 24.8-27.3 cm

MODEL: 2547



ANKLE SUPPORT WITH WRAPAROUND STRAP

HPCPS: A4466 ●● MEDIUM

- Fabricated from neoprene and other lightweight materials to provide medial-lateral stability along with a full range of dorsi- and plantar-flexion in the foot and ankle
- Wide elastic straps wrap around joint in figure-8 style

INDICATIONS: Ankle sprains, chronic ankle instability

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2547-S	S	6.75-7.75 in / 17.1-19.7 cm
2547-M	M	7.75-8.75 in / 19.7-22.2 cm
2547-L	L	8.75-9.75 in / 22.2-24.8 cm
2547-XL	XL	9.75-10.75 in / 24.8-27.3 cm

MODEL: 2560



ANKLE SUPPORT WITH SPIRAL STAYS

HPCPS: A4466 ●● MEDIUM

- Helps reduce swelling and minimize pain while allowing full range of motion.
- Bilateral spiral stays move with the ankle and provide medial-lateral support and stability
- Excellent for light sports and exercise activity

INDICATIONS: Acute ankle sprains, chronic ankle instability, post cast removal

MEASURE AROUND ANKLE AT THE SMALLEST POINT

2560-S	S	6.75-7.75 in / 17.1-19.7 cm
2560-M	M	7.75-8.75 in / 19.7-22.2 cm
2560-L	L	8.75-9.75 in / 22.2-24.8 cm
2560-XL	XL	9.75-10.75 in / 24.8-27.3 cm

MODEL: 2095



LOW PROFILE NIGHT SPLINT WITH AIR PAD

HPCPS: L4396 ●●● MAXIMUM

- Holds the foot in a neutral position during sleep, alleviating the waking pain and discomfort associated with plantar-fasciitis and other conditions
- Lightweight, rigid shell with soft padding and encircling straps for wearing comfort
- Can be worn on either right or left foot

INDICATIONS: Plantar fasciitis, Achilles tendonitis, post-static pain, drop foot

UNIVERSAL FIT

2095-M	M	FOR PETITES, ADOLESCENTS
2095-L	L	FOR AVERAGE TO LARGE ADULTS

MODEL: 2096



SOFT TOP POST-OP SHOE

HPCPS: L3260

- Semi-rigid sole provides superior protection and control
- Soft upper allows the foot to breathe and provides enhanced wearing comfort
- Adjustable straps accommodate extra bulky dressings
- Can be worn on either right or left foot

INDICATIONS: Post foot surgery, post metatarsal injury, soft tissue procedures

SHOE SIZE

2096/M-S	M/S	MEN'S 4 - 6
2096/M-M	M/M	MEN'S 6.5 - 8.5
2096/M-L	M/L	MEN'S 9 - 11
2096/M-XL	M/XL	MEN'S 11.5 - 14
2096/W-S	W/S	WOMEN'S 5 - 6
2096/W-M	W/M	WOMEN'S 7 - 8
2096/W-L	W/L	WOMEN'S 9 - 11

MODEL: 1793



PNEUMATIC WALKER, 17" UPRIGHT

HPCPS: L4360 ●●● MAXIMUM

- Uprights molded onto footplate to form a single piece
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Low and wide rocker bottom helps promote a natural, stable gait
- Accommodates the swelling patterns that occur throughout the rehabilitation process
- Pneumatic bladder is easily adjusted to achieve desired compression and custom fit

INDICATIONS: Soft tissue injuries, grade 2 and 3 sprains, stable fractures, post-operative use, trauma and rehab

SHOE SIZE

		MEN'S	WOMEN'S
1793-S	S	4.5 - 7	6 - 8
1793-M	M	7.5 - 10.5	8.5 - 11.5
1793-L	L	10.5 - 12.5	11.5 - 13.5
1793-XL	XL	13 - 14	14 - 15



ANKLE AND FOOT SUPPORTS



MODEL: 1791



HIGH TOP WALKER, 17" UPRIGHT
HCPCS: L4386 ●●● MAXIMUM

- Uprights molded onto footplate to form a single piece
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Low and wide rocker bottom helps promote a natural, stable gait

INDICATIONS: Soft tissue injuries, grade 2 and 3 sprains, stable fractures, post-operative use, trauma and rehab

SHOE SIZE

		MEN'S	WOMEN'S
1791-S	S	4.5 – 7	6 – 8
1791-M	M	7.5 – 10.5	8.5 – 11.5
1791-L	L	10.5 – 12.5	11.5 – 13.5
1791-XL	XL	13 – 14	14 – 15

MODEL: 1792



LOW TOP WALKER, 11" UPRIGHT
HCPCS: L4386 ●●● MAXIMUM

- Uprights molded onto footplate to form a single piece
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Low and wide rocker bottom helps promote a natural, stable gait
- Lower struts and a low heel height

INDICATIONS: Metatarsal fracture, acute ankle sprains, stress and stable fractures of the foot, severe ankle sprain

SHOE SIZE

		MEN'S	WOMEN'S
1792-XS	XS	2 – 4	3-5
1792-S	S	4.5 – 7	6 – 8
1792-M	M	7.5 – 10.5	8.5 – 11.5
1792-L	L	10.5 – 12.5	11.5 – 13.5
1792-XL	XL	13 – 14	14 – 15

MODEL: CMS-001-17



VALULINE HIGH TOP WALKER, 17" UPRIGHT
HCPCS: L4386 ●●● MAXIMUM

- Uprights riveted onto footplate
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Rocker bottom helps promote a natural, stable gait

INDICATIONS: Soft tissue injuries, grade 2 and 3 sprains, stable fractures, post-operative use, trauma and rehab

SHOE SIZE

		MEN'S	WOMEN'S
CMS-001-17-S	S	4.5 – 7	6 – 8
CMS-001-17-M	M	7.5 – 10.5	8.5 – 11.5
CMS-001-17-M	L	10.5 – 12.5	11.5 – 13.5

MODEL: CMS-001-11



VALULINE LOW TOP WALKER, 11" UPRIGHT
HCPCS: L4386 ●●● MAXIMUM

- Uprights riveted onto footplate
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Rocker bottom helps promote a natural, stable gait
- Lower struts and a low heel height

INDICATIONS: Metatarsal fracture, acute ankle sprains, stress and stable fractures of the foot, severe ankle sprain

SHOE SIZE

		MEN'S	WOMEN'S
CMS-001-11-S	S	4.5 – 7	6 – 8
CMS-001-11-M	M	7.5 – 10.5	8.5 – 11.5
CMS-001-11-M	L	10.5 – 12.5	11.5 – 13.5

MODEL: CMS-003-17



VALULINE PNEUMATIC WALKER, 17" UPRIGHT
HCPCS: L4360 ●●● MAXIMUM

- Uprights riveted onto footplate
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Rocker bottom helps promote a natural, stable gait
- Pneumatic system provides total contact support, increasing stability while also decreasing pain and swelling.

INDICATIONS: Soft tissue injuries, grade 2 and 3 sprains, stable fractures, post-operative use, trauma and rehab

SHOE SIZE

		MEN'S	WOMEN'S
CMS-003-17-S	S	4.5 – 7	6 – 8
CMS-003-17-M	M	7.5 – 10.5	8.5 – 11.5
CMS-003-17-M	L	10.5 – 12.5	11.5 – 13.5

MODEL: CMS-003-11



VALULINE PNEUMATIC WALKER, 11" UPRIGHT
HCPCS: L4360 ●●● MAXIMUM

- Uprights riveted onto footplate
- Shock absorbing sole reduces impact from heel strike to improve comfort during ambulation
- Rocker bottom helps promote a natural, stable gait
- Pneumatic system provides total contact support, increasing stability while also decreasing pain and swelling.
- Lower struts and a low heel height

INDICATIONS: Metatarsal fracture, acute ankle sprains, stress and stable fractures of the foot, severe ankle sprain

SHOE SIZE

		MEN'S	WOMEN'S
CMS-003-11-S	S	4.5 – 7	6 – 8
CMS-003-11-M	M	7.5 – 10.5	8.5 – 11.5
CMS-003-11-M	L	10.5 – 12.5	11.5 – 13.5





ELBOW STRAPS AND SUPPORTS

MODEL: 2089



TENNIS ELBOW STRAP WITH GEL PAD

HCPCS: A4465 •• MEDIUM

- Provides precisely targeted support and compression
- Gel pad is easily detachable from the brush lined encircling strap, and can be cooled or heated for therapy

INDICATIONS: Medial and/or lateral epicondylitis (tennis or golfer's elbow)

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS AND ADOLESCENTS

2089

UNIVERSAL

MODEL: 2421



BANDIT® THERAPEUTIC FOREARM BAND

HCPCS: A4465 •• MEDIUM

- Features two molded plates curved to the shape of the forearm, with soft neoprene cushioning inside
- Can be worn for extended periods without complete comfort without affecting circulation or hampering regular movements

INDICATIONS: Sprains, strains, tendonitis, lateral epicondylitis, carpal tunnel syndrome

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS AND ADOLESCENTS

2421

UNIVERSAL

MODEL: 2423



TENNIS ELBOW STRAP WITH LIGHTWEIGHT AIR PAD

HCPCS: A9270 •• MEDIUM

- Provides precisely targeted support and compression
- Air pad provides focused pneumatic compression that cushions around elbow anatomy

INDICATIONS: Medial and/or lateral epicondylitis (tennis or golfer's elbow)

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS AND ADOLESCENTS

2423

UNIVERSAL

MODEL: 2419



PULLOVER ELASTIC ELBOW SUPPORT

HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the elbow along with all-day wearing comfort
- Cotton-rayon blend fabric helps improve circulation to promote healing

INDICATIONS: Compressive support for weak or injured elbows

MEASURE AROUND THE BEND OF THE ELBOW

2419-S	S	8.5-10 in / 21.6-25.4 cm
2419-M	M	10-11.5 in / 25.4-29.2 cm
2419-L	L	11.5-13 in / 29.2-33 cm
2419-XL	XL	13-14.5 in / 33-36.8 cm

MODEL: 2427



ELBOW SUPPORT WITH VISCOELASTIC INSERT

HCPCS: A4466 •• MEDIUM

- Anatomically sculpted insert pad provides precise compression to support weak or injured elbows
- Controlled, four-way stretch elastic offers effective compressive support to the entire knee area

INDICATIONS: Tendonitis, carpal tunnel syndrome, lateral epicondylitis (tennis elbow), repetitive stress injuries, cumulative trauma disorders

MEASURE AROUND THE BEND OF THE ELBOW

2427-S	S	8.5-10 in / 21.6-25.4 cm
2427-M	M	10-11.5 in / 25.4-29.2 cm
2427-L	L	11.5-13 in / 29.2-33 cm
2427-XL	XL	13-14.5 in / 33-36.8 cm

MODEL: 2460



CRADLE ARM SLING

HCPCS: A4565 •• MEDIUM

- Provides well balanced support to the arm during recuperation from fracture, sprain or surgery, evenly distributing its weight
- Reversible for right or left arm, with wide range of adjustment to easily fits most adults
- Thumb loop supports the hand, maintains the wrist in a neutral position for wearing comfort

INDICATIONS: Cast support, injury to the arm, wrist or hand

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2460

UNIVERSAL



WRIST SUPPORTS



MODEL: 2072



Specify for RIGHT or LEFT hand

THUMB IMMOBILIZER

HPCPS: L3923 •• MEDIUM

- Helps support and limit motion of the thumb
- Allows for full finger dexterity
- Easy fastening and adjustment

INDICATIONS: Soft tissue injuries, ligament strains, Gamekeeper's thumb, osteoarthritis, degenerative joint disease

MEASURE AROUND THE BREAK OF THE WRIST

2072-XS	XS	4.5-5.5 in / 11.4-14 cm
2072-S	S	5.5-6.5 in / 14-16.5 cm
2072-M	M	6.5-7.5 in / 16.5-19 cm
2072-L	L	7.5-8.5 in / 19-21.6 cm
2072-XL	XL	8.5-9.5 in / 21.6-24.1 cm

MODEL: 2081



Specify for RIGHT or LEFT wrist

LOW PROFILE WRIST BRACE

HPCPS: L3908 •• MEDIUM

- Lightweight, low profile exoskeleton provides support with less bulk
- Three dimensionally molded interior creates a padded surface that enhances comfort, relieves predictable pressure points
- Easy hook-and-eye fastening mechanism is adjustable for comfort

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST AND PALM

		WRIST	PALM
2081-S	S	5.25-6.5 in / 13.3-16.5 cm	6.75-8 in / 17.1-20.3 cm
2081-M	M	6.5-7.75 in / 16.5-19.7 cm	8-9.25 in / 20.3-23.5 cm
2081-L	L	7.75-10 in / 19.7-25.4 cm	9.25-10.5 in / 23.5-26.7 cm

MODEL: 2082



Specify for RIGHT or LEFT wrist

6" WRIST SPLINT

HPCPS: L3908 ••• MAXIMUM

- For petite to average females and adolescents
- Well-padded lightweight glove with separate cinch fasteners provides customized support and compression
- Low contour at palmar crease allows for full finger dexterity and grip strength
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2082-XS	XS	4.5-5.5 in / 11.4-14 cm
2082-S	S	5.5-6.5 in / 14-16.5 cm
2082-M	M	6.5-7.5 in / 16.5-19 cm
2082-L	L	7.5-8.5 in / 19-21.6 cm
2082-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2083



Specify for RIGHT or LEFT wrist

8" WRIST SPLINT

HPCPS: L3908 ••• MAXIMUM

- For average to larger adults
- Well-padded lightweight glove with separate cinch fasteners provides customized support and compression
- Low contour at palmar crease allows for full finger dexterity and grip strength
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2083-XS	XS	4.5-5.5 in / 11.4-14 cm
2083-S	S	5.5-6.5 in / 14-16.5 cm
2083-M	M	6.5-7.5 in / 16.5-19 cm
2083-L	L	7.5-8.5 in / 19-21.6 cm
2083-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2086



Specify for RIGHT or LEFT wrist

6" WRIST-THUMB SPLINT

HPCPS: L3908 ••• MAXIMUM

- For petite to average females and adolescents
- Well-padded lightweight glove with separate cinch fasteners provides customized support and compression
- Radial splint anchored by proximal circumference strap provides for optimum immobilization of the thumb
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: DeQuervain's tendonitis, strains and sprains of the MCP joint, basal joint arthritis, advanced carpal tunnel syndrome, cumulative trauma disorders of the thumb and wrist

MEASURE AROUND THE BREAK OF THE WRIST

2086-S	S	5.5-6.5 in / 14-16.5 cm
2086-M	M	6.5-7.5 in / 16.5-19 cm
2086-L	L	7.5-8.5 in / 19-21.6 cm
2086-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2087



Specify for RIGHT or LEFT wrist

8" WRIST-THUMB SPLINT

HPCPS: L3908 ••• MAXIMUM

- For average to larger adults
- Well-padded lightweight glove with separate cinch fasteners provides customized support and compression
- Radial splint anchored by proximal circumference strap provides for optimum immobilization of the thumb
- Removable palmar splint and flexible stays for maximum support and immobilization

INDICATIONS: DeQuervain's tendonitis, strains and sprains of the MCP joint, basal joint arthritis, advanced carpal tunnel syndrome, cumulative trauma disorders of the thumb and wrist

MEASURE AROUND THE BREAK OF THE WRIST

2087-XS	XS	4.5-5.5 in / 11.4-14 cm
2087-S	S	5.5-6.5 in / 14-16.5 cm
2087-M	M	6.5-7.5 in / 16.5-19 cm
2087-L	L	7.5-8.5 in / 19-21.6 cm
2087-XL	XL	8.5-10 in / 21.6-25.4 cm





WRIST SUPPORTS

MODEL: 2351



REVERSIBLE COCK-UP WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Reversible to fit either right or left wrist
- Lined elastic glove provides support and compression, wicking away moisture for long term wearing comfort
- Opens flat for easy application over injured or tender wrists

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2351-XS	XS	4.5-5.5 in / 11.4-14 cm
2351-S	S	5.5-6.5 in / 14-16.5 cm
2351-M	M	6.5-7.5 in / 16.5-19 cm
2351-L	L	7.5-8.5 in / 19-21.6 cm
2351-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 2360



SUEDE FINISH COCK-UP WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Durable suede finish glove is perforated and flannel lined for comfortable, long term support
- Lacing closure tab provides uniform compression and support
- Preshaped splint is removable for additional contouring if desired

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2360-XS	XS	3.5-4.5 in / 8.9-11.4 cm
2360-S	S	4.5-5.5 in / 11.4-14 cm
2360-M	M	5.5-6.5 in / 14-16.5 cm
2360-L	L	6.5-7.5 in / 16.5-19 cm
2360-XL	XL	7.5-8.5 in / 19-21.6 cm

Specify for RIGHT or LEFT wrist

MODEL: 2362



CANVAS COCK-UP WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Durable canvas glove is flannel lined for comfortable, full-length support
- Criss-cross lacing provides uniform compression and support
- Preshaped splint is removable for additional contouring if desired

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2362-2XS	2XS	3.5-4.5 in / 8.9-11.4 cm
2362-XS	XS	4.5-5.5 in / 11.4-14 cm
2362-S	S	5.5-6.5 in / 14-16.5 cm
2362-M	M	6.5-7.5 in / 16.5-19 cm
2362-L	L	7.5-8.5 in / 19-21.6 cm
2362-XL	XL	8.5-9.5 in / 21.6-24.1 cm

Specify for RIGHT or LEFT wrist

MODEL: 2364



PROFESSIONALS CHOICE COCK-UP WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Available in durable white [no. 2364] or gray canvas [no. 2364E]
- Canvas glove is flannel lined for comfortable, full-length support
- Criss-cross lacing provides uniform compression and support
- Preshaped splint is removable for additional contouring if desired

INDICATIONS: Carpal tunnel syndrome, mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

2364/E-2XS	2XS	3.5-4.5 in / 8.9-11.4 cm
2364/E-XS	XS	4.5-5.5 in / 11.4-14 cm
2364/E-S	S	5.5-6.5 in / 14-16.5 cm
2364/E-M	M	6.5-7.5 in / 16.5-19 cm
2364/E-L	L	7.5-8.5 in / 19-21.6 cm
2364/E-XL	XL	8.5-9.5 in / 21.6-24.1 cm

Specify for RIGHT or LEFT wrist

MODEL: 2389



OCCUPATIONAL WRIST SUPPORT HCPCS: A4466 •• MEDIUM

- Excellent choice for workplace use to prevent and/or treat carpal tunnel syndrome
- Lightweight elastic glove and flexible palmar stays provide support to the wrist while permitting free motion and full finger dexterity

INDICATIONS: For wear while performing highly repetitive manual tasks [data entry, etc.]

MEASURE AROUND THE BREAK OF THE WRIST

2389-2XS	2XS	3.5-4.5 in / 8.9-11.4 cm
2389-XS	XS	4.5-5.5 in / 11.4-14 cm
2389-S	S	5.5-6.5 in / 14-16.5 cm
2389-M	M	6.5-7.5 in / 16.5-19 cm
2389-L	L	7.5-8.5 in / 19-21.6 cm
2389-XL	XL	8.5-9.5 in / 21.6-24.1 cm

Specify for RIGHT or LEFT wrist

MODEL: 2418



PULLOVER ELASTIC WRIST SUPPORT HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the wrist along with all-day wearing comfort
- Cotton-rayon blend fabric helps improve circulation to promote healing

INDICATIONS: Compressive support for weak or injured wrists

MEASURE AROUND THE BREAK OF THE WRIST

2418-S	S	5.5-6.5 in / 14-16.5 cm
2418-M	M	6.5-7.5 in / 16.5-19 cm
2418-L	L	7.5-8.5 in / 19-21.6 cm
2418-XL	XL	8.5-9.5 in / 21.6-24.1 cm



MATERNITY AND BACK SUPPORTS



MODEL: 2786



MATERNITY SUPPORT

HPCPS: A9270 •• MEDIUM

- Provides effective uplift support over the extended abdomen, stabilizing the ligaments, tendons and muscles of the lower back, and offering relief from strain and discomfort
- Made from finest quality beige hinged elastic that molds comfortably around surface anatomy.
- Flexible boning in back panel provides additional support

INDICATIONS: Relief from lower back discomfort during pregnancy

SIZED ACCORDING TO DRESS SIZE
PRIOR TO PREGNANCY

2786-S	S	DRESS SIZE 2 - 6
2786-M	M	DRESS SIZE 6 - 14
2786-L	L	DRESS SIZE 16 - 20

MODEL: 2892



PETITE LUMBOSACRAL SUPPORT

HPCPS: A4466 •• MEDIUM

- Effective stabilizes soft tissue support structures of the lower back, providing relief from pain and discomfort
- Made from finest quality beige hinged elastic that molds comfortably around surface anatomy
- Flexible boning in back panel provides additional support
- 6" deep in front, 9" deep in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2892-S	S	29-32 in / 73.6-81.3 cm
2892-M	M	33-36 in / 81.3-91.4 cm
2892-L	L	37-40 in / 94-101.6 cm
2892-XL	XL	41-44 in / 104.1-111.8 cm
2892-U	UNIVERSAL	29-44 in / 73.6-111.8 cm

MODEL: 2893



LUMBOSACRAL SUPPORT WITH ABDOMINAL UPLIFT

HPCPS: A4466 ••• MAXIMUM

- Provides firm, form-fitting support for relief of low back pain and discomfort
- Strong elastic abdominal belt provides additional uplift support and is adjustable for varying degree of compression
- Two rigid metal stays are shapeable for extra firm support
- 7" deep in front, 12" in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2893-S	S	29-32 in / 73.6-81.3 cm
2893-M	M	33-36 in / 81.3-91.4 cm
2893-L	L	37-40 in / 94-101.6 cm
2893-XL	XL	41-44 in / 104.1-111.8 cm
2893-U	UNIVERSAL	29-44 in / 73.6-111.8 cm

MODEL: 2884



LIGHTWEIGHT LUMBOSACRAL SUPPORT

HPCPS: A4466 ••• MAXIMUM

- Strong mesh elastic material conforms to surface anatomy for a cool, comfortable custom fit
- Optional metal stays can be shaped to provide extra-firm support to the spine
- Wide elastic adjustment panels
- 7" deep in front, 11" deep in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2884-S	S	26-32 in / 66-81.3 cm
2884-M	M	32-38 in / 81.3-96.5 cm
2884-L	L	38-44 in / 96.5-111.8 cm
2884-XL	XL	44-50 in / 111.8-127 cm
2884-U	UNIVERSAL	26-44 in / 66-111.8 cm

MODEL: 2890



LUMBOSACRAL SUPPORT WITH ABDOMINAL UPLIFT

HPCPS: A4466 ••• MAXIMUM

- Provides firm support and stabilization to aid in relief of low back pain and discomfort
- Strong elastic abdominal belt provides additional uplift support and is adjustable for varying degree of compression
- Two rigid metal stays are shapeable for extra firm support
- 7" deep in front, 11" in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2890-S	S	29-32 in / 73.6-81.3 cm
2890-M	M	33-36 in / 81.3-91.4 cm
2890-L	L	37-40 in / 94-101.6 cm
2890-XL	XL	41-44 in / 104.1-111.8 cm
2890-2XL	2XL	41-44 in / 104.1-111.8 cm
2890-U	UNIVERSAL	29-44 in / 73.6-111.8 cm

MODEL: 2891



PETITE LIGHTWEIGHT LUMBOSACRAL SUPPORT

HPCPS: A4466 •• MEDIUM

- Designed for the very petite female figure
- Lightweight ribbed elastic provides form-fitting support to relieve lower back discomfort and pain
- Flexible stays provide added support and stability
- 7" deep back panel tapers to 4.5" in front

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2891-S	S	29-32 in / 73.6-81.3 cm
2891-M	M	33-36 in / 81.3-91.4 cm
2891-L	L	37-40 in / 94-101.6 cm
2891-XL	XL	41-44 in / 104.1-111.8 cm
2891-U	UNIVERSAL	29-44 in / 73.6-111.8 cm





BACK SUPPORTS

MODEL: 2886



HEAVY DUTY SACRO BRACE WITH THERMOPAD

HPCPS: LA4466 •• MEDIUM

- Unique construction provides optimum support and comfort to the sacrum and the lower spine
- Adjustable tapered design conforms to effectively support virtually every adult figure type, male or female
- Removable ThermoPad provides added warmth for aching muscles, joints, tendons and discs that so often accompany low back pain
- 6" deep in front, 9" deep in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2886 UNIVERSAL 28-50 in

MODEL: 2895



LUMBOTeK LUMBOSACRAL SUPPORT

HPCPS: L0642 ••• MAXIMUM

- Extends from lower sacrum to L-1 vertebra
- Made from high quality elastic and cotton blend materials
- Rigid anterior and posterior panels provide additional support to the abdomen and relieve pressure over the lower back
- Varying degree of additional support can be obtained by adjusting the elastic support pulls

INDICATIONS: Lumbar muscle or ligament sprain or strain, chronic low back pain, herniated disc or disc disease, wherever unloading of disc pressure is desired, compression fracture, spinal stabilization, post laminectomy, spondylolysis, spondylolisthesis, posture support

MEASURE AROUND THE
WIDEST PART OF THE HIPS

		WAIST	DEPTH	COLOR
2895BL09-S	S	26 – 32"	9"	BLACK
2895BL09-M	M	30 – 37"	9"	BLACK
2895BL09-L	L	34 – 41"	9"	BLACK
2895BL09-XL	XL	38 – 47"	9"	BLACK
2895BL09-2XL	2 XL	44 – 54"	9"	BLACK
2895BL09-S3XL	3 XL	50 – 62"	9"	BLACK
2895BL11-S	S	26 – 32"	11"	BLACK
2895BL11-M	M	30 – 37"	11"	BLACK
2895BL11-L	L	34 – 41"	11"	BLACK
2895BL11-XL	XL	38 – 47"	11"	BLACK
2895BL11-2XL	2 XL	44 – 54"	11"	BLACK
2895BL11-3XL	3 XL	50 – 62"	11"	BLACK
2895WH09-S	S	26 – 32"	9"	WHITE
2895WH09-M	M	30 – 37"	9"	WHITE
2895WH09-L	L	34 – 41"	9"	WHITE
2895WH09-XL	XL	38 – 47"	9"	WHITE
2895WH09-2XL	2 XL	44 – 54"	9"	WHITE
2895WH09-3XL	3 XL	50 – 62"	9"	WHITE

MODEL: 2894



Showing the interior of an unfastened 2894 TruTEK LS orthosis

TRUTEK LUMBOSACRAL SUPPORT

HPCPS: L0648 ••• MAXIMUM

- Extends from lower sacrum to T-9 vertebra
- Made from high quality nylon and cotton blend materials
- Rigid anterior and posterior panels provide additional support to the abdomen and relieve pressure over the lower back
- Varying degree of additional support can be obtained by adjusting the elastic support pulls
- 14" deep in back, 8" deep in front

INDICATIONS: Lumbar muscle or ligament sprain or strain, chronic low back pain, herniated disc or disc disease, wherever unloading of disc pressure is desired, compression fracture, spinal stabilization, post laminectomy, spondylolysis, spondylolisthesis, posture support

MEASURE AROUND THE
WIDEST PART OF THE HIPS

2894-S	S	26-32 in / 66-81.3 cm
2894-M	M	30-37 in / 76.2-94 cm
2894-L	L	34-41 in / 86.4-104.1 cm
2894-XL	XL	38-47 in / 96.5-119.4 cm
2894-2XL	2XL	44-54 in / 111.8-137.2 cm
2894-3XL	3XL	50-62 in / 127-157.5 cm



RIB SUPPORTS



MODEL: 2459



RIB BELT FOR MEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Wraparound elastic web material with adjustable contact closure fastening.
- 6" deep all around

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE THE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2459-UR	UR	29 – 38 in / 73.6-96.5 cm
2459-UL	UL	39 – 48 in / 99-121.9 cm

MODEL: 2659



RIB BELT FOR WOMEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Wraparound elastic web material with adjustable contact closure fastening.
- Cut low in front to accommodate the bust line
- 4" in front, 6" in back

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE THE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2659-UR	UR	29 – 38 in / 73.6-96.5 cm
2659-UL	UL	39 – 48 in / 99-121.9 cm

ABDOMINAL AND HERNIA SUPPORTS



MODEL: 2509



THREE PANEL ELASTIC BINDER FOR MEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Hinged three-panel elastic molds around surface anatomy for a comfortable, custom fit
- Can also be used as a rib fracture support
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2509-S	S	29-32 in / 73.6-81.3 cm
2509-M	M	33-36 in / 81.3-91.4 cm
2509-L	L	37-40 in / 94-101.6 cm
2509-XL	XL	41-44 in / 104.1-111.8 cm
2509-2XL	2XL	45-48 in / 114.3-121.9 cm

MODEL: 2689



THREE PANEL ELASTIC BINDER FOR WOMEN

HPCPS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Hinged three-panel elastic molds around surface anatomy for a comfortable, custom fit
- Contoured to fit the female form
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2689-S	S	19-22 in / 48.3-55.9 cm
2689-M	M	23-26 in / 58.4-66 cm
2689-L	L	27-30 in / 68.6-76.2 cm
2689-XL	XL	31-34 in / 78.7-86.4 cm
2689-2XL	2XL	35-38 in / 88.9-96.5 cm

MODEL: 2518



MULTIPLE USE ELASTIC ABDOMINAL BINDER

HPCPS: A4466 • MINIMUM

- Lightweight full length binder has a variety of applications, most commonly to protect post-surgical dressings
- Effective in supporting abdominal contents when there has been a loss of muscle tone
- Flexible plastic support stays prevent rolling
- 10" deep all around

INDICATIONS: Umbilical hernia, weak or pendulous abdomen, abdominal contusions and tenderness, minor rib fracture

MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2518-S	S	29-32 in / 73.6-81.3 cm
2518-M	M	33-36 in / 81.3-91.4 cm
2518-L	L	37-40 in / 94-101.6 cm
2518-XL	XL	41-44 in / 104.1-111.8 cm
2518-2XL	2XL	45-48 in / 114.3-121.9 cm

MODEL: 2956



LIGHTWEIGHT BEIGE HERNIA SUPPORT

HPCPS: L8300 (SINGLE) •• MEDIUM
L8310 (BILATERAL)

- Plush lightweight material is contact closure compatible
- Elastic waistband in back molds around surface anatomy for a comfortable, custom fit
- Removable foam pads provide gentle, localized support at the site of the hernia
- Soft elastic understraps can be tightened or loosened to increase or lessen support
- Completely free of any metal component

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST PART OF THE HIPS

2956-S	S	30-35 in / 76-89 cm
2956-M	M	35-40 in / 89-102 cm
2956-L	L	41-46 in / 104-117 cm





HERNIA SUPPORTS

MODEL: 2936 SINGLE MODEL: 2939 BILATERAL



HERNIA SUPPORT (SCROTAL ELASTIC TRUSS)

HCPCS: L8300 (2936) •• MEDIUM
L8310 (2939)

- Relieves discomfort while treating inguinal or scrotal hernias
- Design of the truss allows the specially shaped pad to self-align and maintain its proper position, thus applying support as a constant, unified force.

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST PART OF THE HIPS

2936/L-S	S	LEFT	30-32 in / 76.2-81.3 cm
2936/L-M	M	LEFT	34-36 in / 86.4-91.4 cm
2936/L-L	L	LEFT	38-40 in / 96.5-101.6 cm
2936/L-XL	XL	LEFT	42-44 in / 106.7-111.8 cm
2936/L-2XL	2 XL	LEFT	45-48 in / 114.3-121.9 cm
2936/R-S	S	RIGHT	30-32 in / 76.2-81.3 cm
2936/R-M	M	RIGHT	34-36 in / 86.4-91.4 cm
2936/R-L	L	RIGHT	38-40 in / 96.5-101.6 cm
2936/R-XL	XL	RIGHT	42-44 in / 106.7-111.8 cm
2936/R-2XL	2 XL	RIGHT	45-48 in / 114.3-121.9 cm
2939-S	S	BILATERAL	30-32 in / 76.2-81.3 cm
2939-M	M	BILATERAL	34-36 in / 86.4-91.4 cm
2939-L	L	BILATERAL	38-40 in / 96.5-101.6 cm
2939-XL	XL	BILATERAL	42-44 in / 106.7-111.8 cm
2939-2XL	2 XL	BILATERAL	45-48 in / 114.3-121.9 cm

MODEL: 2958



Specify for hernia on wearer's RIGHT or LEFT side

LIGHTWEIGHT SINGLE HERNIA SUPPORT

HCPCS: L8300 •• MEDIUM

- Holds a single inguinal or scrotal hernia comfortably in reduction
- Can also be used in special situations after a hernia repair to prevent recurrence
- Low profile asymmetric design allows it to be worn under clothing without detection
- 4" web elastic encircling belt molds around body contours for a comfortable, custom fit

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST PART OF THE HIPS

2958-S	S	30-32 in / 76.2-81.3 cm
2958-M	M	34-36 in / 86.4-91.4 cm
2958-L	L	38-40 in / 96.5-101.6 cm
2958-XL	XL	42-44 in / 106.7-111.8 cm

MODEL: 2958D



LIGHTWEIGHT DOUBLE HERNIA SUPPORT

HCPCS: L8310 •• MEDIUM

- Relieves discomfort while treating inguinal or scrotal hernias
- Can also be used after hernia repair to prevent recurrence
- Removable foam pads provide localized support at hernia site
- Detachable understraps can be tightened or loosened to increase or lessen support

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST PART OF THE HIPS

2958D-S	S	30-32 in / 76.2-81.3 cm
2958D-M	M	34-36 in / 86.4-91.4 cm
2958D-L	L	38-40 in / 96.5-101.6 cm
2958D-XL	XL	42-44 in / 106.7-111.8 cm



CERVICAL AND CLAVICLE SUPPORTS

MODEL: 2452



POSTURE SUPPORT

HCPCS: A4466 • MINIMUM

- Serves as a gentle reminder to hold the shoulders in a natural, correct posture, especially while performing highly repetitive tasks like data entry
- Also helps prevent strain and fatigue
- Made from finest quality beige hinged elastic web: 3" wide shoulder straps, 6" back panel

INDICATIONS: Slumping shoulders, poor posture

MEASURE WAIST CIRCUMFERENCE

2452-XS	XS	25-28.5 in / 63-72 cm
2452-S	S	28.5-31.5 in / 72-80 cm
2452-M	M	31.5-35.5 in / 80-90 cm
2452-L	L	35.5-39.5 in / 90-100 cm
2452-XL	XL	39.5-43.5 in / 100-110 cm

MODEL: 2455



POSTURE CONTROL SHOULDER BRACE

HCPCS: L3650 •• MEDIUM

- Serves as a gentle reminder to hold the shoulders in a natural, correct posture, especially while performing highly repetitive tasks like data entry
- Also helps prevent strain and fatigue
- Comfortable white canvas back panel, non-elastic web shoulder straps

INDICATIONS: Slumping shoulders, poor posture

MEASURE AROUND CHEST AT THE BASE OF THE BREAST BONE

2455-XS	XS	UP TO 26 in / UP TO 66 cm
2455-S	S	27-32 in / 68.9-81.3 cm
2455-M	M	33-38 in / 83.8-96.5 cm
2455-L	L	39-44 in / 99-111.7 cm
2455-XL	XL	OVER 45 in / OVER 114.3 cm



CERVICAL AND CLAVICLE SUPPORTS



MODEL: 2453



CLAVICLE STRAP

HCPCS: L3650 •• MEDIUM

- Holds the shoulder region in the correct position to help promote healing following a clavicle fracture or separation
- Figure 8 design for easy application
- Varying degrees of support can be obtained by adjusting the tension of the shoulder straps

INDICATIONS: Clavicle fracture or separation, poor posture

MEASURE CHEST CIRCUMFERENCE AT THE LEVEL OF THE UNDERARMS

2453-XS	XS	20-24 in / 51-61 cm
2453-S	S	24-30 in / 61-76 cm
2453-M	M	30-36 in / 76-91 cm
2453-L	L	36-42 in / 91-107 cm
2453-XL	XL	42-48 in / 107-122 cm

MODEL: 2394



SOFT FOAM CERVICAL COLLAR

HCPCS: L0120 • MINIMUM

- Medium density foam covered in stockinette
- Especially convenient for wear while sleeping
- Narrow, average and wide depths

INDICATIONS: Relief from muscle tension, treatment of minor neck pain or injuries, whiplash, as a reminder against making sudden and painful neck movements

MEASURE NECK CIRCUMFERENCE

	NECK	DEPTH*
2394/N-XS	11.5-13.5 in / 29.2-34.3 cm	2.5"
2394/N-S	13.5-15.5 in / 34.3-39.4 cm	2.5"
2394/N-M	15.5-17.5 in / 39.4-44.4 cm	2.5"
2394/N-L	17.5-19.5 in / 44.4-49.5 cm	2.5"
2394/N-XL	19.5-21.5 in / 49.5-54.6 cm	2.5"
2394/N-U	13.5-19.5 in / 34.3-49.5 cm	2.5"
2394/A-XS	11.5-13.5 in / 29.2-34.3 cm	3"
2394/A-S	13.5-15.5 in / 34.3-39.4 cm	3"
2394/A-M	15.5-17.5 in / 39.4-44.4 cm	3"
2394/A-L	17.5-19.5 in / 44.4-49.5 cm	3"
2394/A-XL	19.5-21.5 in / 49.5-54.6 cm	3"
2394/A-U	13.5-19.5 in / 34.3-49.5 cm	3"
2394/W-XS	11.5-13.5 in / 29.2-34.3 cm	3.5"
2394/W-S	13.5-15.5 in / 34.3-39.4 cm	3.5"
2394/W-M	15.5-17.5 in / 39.4-44.4 cm	3.5"
2394/W-L	17.5-19.5 in / 44.4-49.5 cm	3.5"
2394/W-XL	19.5-21.5 in / 49.5-54.6 cm	3.5"
2394/W-U	13.5-19.5 in / 34.3-49.5 cm	3.5"

*For front depth, take straight line measurements from sternal notch to base of mandible

MODEL: 2501



OVER-THE-DOOR CERVICAL TRACTION KIT

HCPCS: E0860 ••• MAXIMUM

- Relieves pressure on muscle, nerve tissue and minor fractures of the cervical spine
- Is easily set up for home use
- Includes head halter, spreader bar, door bracket, door hanger and two pulleys, section of rope and water weight bag
- Weight bag is calibrated from 2 lbs to 20 lbs; easy to fill with ordinary tap water

INDICATIONS: Nerve pain from herniated discs, vertebral joint pain, muscle spasm, loss of spinal motion
CONTRAINDICATIONS: Acute cervical injury, osteomyelitis, spinal instability and hypermobility, tumors, rheumatoid arthritis

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2501 UNIVERSAL

MODEL: 2466



GAIT/TRANSFER BELT

HCPCS: A4565 •• MEDIUM

- 2" wide heavy duty web belt with locking teeth fastening mechanism
- Helps facilitate patient transfers
- Decreases the chance of injury to the attendant while helping a patient move or walk

INDICATIONS: Used for moving a person with balance problems from one place to another; to assist a person with balance problems in walking

AVAILABLE IN TWO LENGTHS

2466-48	48 in / 121.9 cm
2466-54	54 in / 137.1 cm

MODEL: 2465



SLING AND SWATHE SHOULDER IMMOBILIZER

HCPCS: L3670 ••• MAXIMUM

- Restricts movement in and around the shoulder by reducing abduction and arm rotation
- Lightweight materials, generous hook and loop closure and adjustment
- Simple to apply, comfortable to wear, and suitable for most adults

INDICATIONS: Shoulder dislocations, separations, minor rotator cuff injuries, strains, post-traumatic, post-surgical conditions

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2465 UNIVERSAL

MODEL: 2464



LIGHTWEIGHT SHOULDER IMMOBILIZER

HCPCS: A4565 •• MEDIUM


- Restricts shoulder and upper arm mobility
- Lightweight, breathable components provide excellent immobilization and wearing comfort
- Removable webbing swathe strap

INDICATIONS: Minor rotator cuff injuries, strains, post-traumatic, post-surgical conditions

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

2464 UNIVERSAL





ARTHRITIC WARMERS

ANGORA ARTHRITIC WARMERS:

- ANGORA WOOL is the lightest of all natural fibers.
- The hollow structure of its fibers makes it ideal for thermal clothing.
- It is comfortable to wear next to the skin, allowing the skin to breathe but absorbing moisture and maintaining a steady temperature at the skin surface.
- ANGORA WOOL have beneficial therapeutic effects for arthritis, carpal tunnel syndrome, back problems and any other pain affecting the joints.

MODEL: 79010



KNEE WARMERS SOLD IN PAIRS

MEASURE AROUND THE BEND OF THE KNEE

79010-S	S	14.5-16 in / 36.8-40.6 cm
79010-M	M	16-18 in / 40.6-45.7 cm
79010-L	L	18-19.5 in / 45.7-49.5 cm
79010-XL	XL	19.5-21 in / 49.5-53.3 cm

MODEL: 79020



LOWER BACK WARMER

MEASURE AROUND THE WAIST

79020-S	S	27-30 in / 68.5-79 cm
79020-M	M	31-35 in / 79-89 cm
79020-L	L	36-40 in / 91.5-101.5 cm
79020-XL	XL	41-45 in / 104-114 cm

MODEL: 79030



SHOULDER WARMER

MEASURE AROUND CHEST

79030-S	S	30-38 in / 76-86 cm
79030-M	M	34-38 in / 86-96.5 cm
79030-L	L	38-42 in / 96.5-107 cm
79030-XL	XL	42-46 in / 107-117 cm

MODEL: 79040



ELBOW / ARM WARMERS SOLD IN PAIRS

MEASURE AROUND THE BEND OF THE ELBOW

79040-S	S	7.75-9.75 in / 19.7-24.7 cm
79040-M	M	9.75-11 in / 24.7-27.9 cm
79040-L	L	11-13.5 in / 27.9-34.3 cm
79040-XL	XL	13.5-15 in / 34.3-38.1 cm

MODEL: 79600



ANGORA SOCKS SOLD IN PAIRS

SHOE SIZE

		MEN'S	WOMEN'S
79600-S	S	N/A	4 - 5
79600-M	M	5 - 7	5.5 - 7.5
79600-L	L	7 - 9.5	8 - 10.5
79600-XL	XL	10 - 11.5	N/A

MODEL: 2088



ARTHRITIS GLOVES

HCPCS: A9270
SOLD IN PAIRS

- Neoprene material retains natural heat to provide soothing warmth and long lasting pain relief
- Interior lining keeps the skin ventilated and comfortably dry

INDICATIONS: Finger and hand pain associated with arthritis, aching wrists, carpal tunnel syndrome, tendonitis

MEASURE AROUND THE HAND
AT THE BASE OF THE KNUCKLES

2088-XS	XS	6.5-7 in / 17.1-17.8 cm
2088-S	S	7-7.75 in / 17.8-19.7 cm
2088-M	M	8-8.75 in / 20.3-22.2 cm
2088-L	L	9-10.5 in / 22.9-26.7 cm
2088-XL	XL	10.75-11.5 in / 27.3-29.2 cm





MODEL: 0141



KNEE SUPPORT WITH OVAL PAD HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over all surface anatomy prominences of the joint
- Extra padding in front for twice the protection against traumatic contact and friction burns over the patellar area

INDICATIONS: Sprains and strains, osteoarthritis

MEASURE AROUND THE BEND OF THE KNEE

0141-S	S	10-12.75 in / 25.4-32.4 cm
0141-M	M	13-15.75 in / 33-40 cm
0141-L	L	16-18.75 in / 40.6-47.6 cm
0141-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0142



KNEE SUPPORT WITH HOR-SHU PAD HCPCS: A4466 • MINIMUM

- Hor-Shu support pad surrounds kneecap and provides compression and stabilization
- Open patella area helps reduce binding around the back of the knee

INDICATIONS: Patellar tendon sprains and strains, post-operative conditions, post-trauma discomfort, inflammation, bursitis

MEASURE AROUND THE BEND OF THE KNEE

0142-S	S	10-12.75 in / 25.4-32.4 cm
0142-M	M	13-15.75 in / 33-40 cm
0142-L	L	16-18.75 in / 40.6-47.6 cm
0142-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0143



KNEE SUPPORT WITH HINGED BARS HCPCS: L1812 ••• MAXIMUM

- Lightweight rigid metal hinged bars provide excellent side to side control and allow for free flexion and extension
- Hor-Shu support pad surrounds kneecap and provides compression and stabilization
- Thigh and calf straps are adjustable to help brace maintain its proper position

INDICATIONS: Collateral ligament damage, patellar tendon sprains and strains, post-operative conditions, post-trauma discomfort

MEASURE AROUND THE BEND OF THE KNEE

0143-S	S	10-12.75 in / 25.4-32.4 cm
0143-M	M	13-15.75 in / 33-40 cm
0143-L	L	16-18.75 in / 40.6-47.6 cm
0143-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0144



KNEE SUPPORT WITH SPIRAL STAYS HCPCS: A4466 •• MEDIUM

- Spiral stays provide moderate medial/lateral stability, yet flex and extend naturally with the joint.
- Hor-Shu support pad surrounds kneecap and provides compression and stabilization
- Thigh and calf straps are adjustable to help brace maintain its proper position

INDICATIONS: Chondromalacia, patellar misalignment and instability, cruciate ligament support, patellar tendonitis

MEASURE AROUND THE BEND OF THE KNEE

0144-S	S	10-12.75 in / 25.4-32.4 cm
0144-M	M	13-15.75 in / 33-40 cm
0144-L	L	16-18.75 in / 40.6-47.6 cm
0144-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0306



NEOPRENE KNEE SUPPORT WITH OPEN PATELLA HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over the surface anatomy of the joint
- Multidirectional stretch materials provide even compression over surface anatomy
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0306RB) or black (0306BL)

INDICATIONS: Sprains and strains, osteoarthritis

MEASURE AROUND THE BEND OF THE KNEE

0306-S	S	12.5-13.75 in / 31.7-34.9 cm
0306-M	M	14-15.25 in / 35.5-38.7 cm
0306-L	L	15.5-16.5 in / 39.3-41.9 cm
0306-XL	XL	16-17.5 in / 40.6-44.4 cm
0306-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
0306-3XL	3XL	19.5-21 in / 49.5-53.3 cm
0306-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 0309



NEOPRENE KNEE SUPPORT WITH STABILIZER PAD HCPCS: A4466 • MINIMUM

- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Multidirectional stretch materials provide even compression over surface anatomy
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0309RB) or black (0309BL)

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

0309-S	S	12.5-13.75 in / 31.7-34.9 cm
0309-M	M	14-15.25 in / 35.5-38.7 cm
0309-L	L	15.5-16.5 in / 39.3-41.9 cm
0309-XL	XL	16-17.5 in / 40.6-44.4 cm
0309-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
0309-3XL	3XL	19.5-21 in / 49.5-53.3 cm
0309-4XL	4XL	21-24 in / 53.3-61 cm





KNEE SUPPORTS

MODEL: 0308



NEOPRENE KNEE STABILIZER WITH SPIRAL STAYS

HCPCS: A4466 •• MEDIUM

- Spiral stays provide moderate medial/lateral stability
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0308RB) or black (0308BL)

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

0308-S	S	12.5-13.75 in / 31.7-34.9 cm
0308-M	M	14-15.25 in / 35.5-38.7 cm
0308-L	L	15.5-16.5 in / 39.3-41.9 cm
0308-XL	XL	16-17.5 in / 40.6-44.4 cm
0308-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
0308-3XL	3XL	19.5-21 in / 49.5-53.3 cm
0308-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 0311



NEOPRENE STABILIZER WRAP WITH SPIRAL STAYS

HCPCS: A4466 •• MEDIUM

- Opens flat for easy application over injured, swollen or tender knees
- Spiral stays provide moderate medial/lateral stability
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0311RB) or black (0311BL)

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

MEASURE AROUND THE BEND OF THE KNEE

0311-XS	XS	12-13 in / 20.5-33 cm
0311-S	S	13-14 in / 33-35.5 cm
0311-M	M	14-15 in / 35.5-38.1 cm
0311-L	L	15-16 in / 38.1-40.6 cm
0311-XL	XL	16-17 in / 40.6-43.2 cm
0311-2XL	2XL	17-18 in / 43.2-45.7 cm
0311-3XL	3XL	18-19 in / 45.7-48.3 cm
0311-4XL	4XL	19-20 in / 48.3-50.8 cm
0311-5XL	5XL	20-22 in / 50.8-55.9 cm

MODEL: 0310



NEOPRENE KNEE STABILIZER WITH HINGED BARS

HCPCS: L1812 ••• MAXIMUM

- Hinged bars provide optimum medial/lateral stability
- Multifunction buttress surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0310RB) or black (0310BL)

INDICATIONS: Medial/lateral knee instability, Chondromalacia, osteoarthritis, compression and support

MEASURE AROUND THE BEND OF THE KNEE

0310-S	S	12.5-13.75 in / 31.7-34.9 cm
0310-M	M	14-15.25 in / 35.5-38.7 cm
0310-L	L	15.5-16.5 in / 39.3-41.9 cm
0310-XL	XL	16-17.5 in / 40.6-44.4 cm
0310-2XL	2XL	17.75-19.25 in / 45.1-48.9 cm
0310-3XL	3XL	19.5-21 in / 49.5-53.3 cm
0310-4XL	4XL	21-24 in / 53.3-61 cm

MODEL: 0312



NEOPRENE STABILIZER WRAP WITH HINGED BARS

HCPCS: L1812 ••• MAXIMUM

- Opens flat for easy application over injured, swollen or tender knees
- Hinged bars provide optimum medial/lateral stability
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Encircling straps are reversible for wear on right or left knee
- Plush lined interior absorbs moisture, wicks away perspiration
- Available in royal (0312RB) or black (0312BL)

INDICATIONS: Medial/lateral knee instability, Chondromalacia, osteoarthritis, compression and support

MEASURE AROUND THE BEND OF THE KNEE

0312-XS	XS	12-13 in / 20.5-33 cm
0312-S	S	13-14 in / 33-35.5 cm
0312-M	M	14-15 in / 35.5-38.1 cm
0312-L	L	15-16 in / 38.1-40.6 cm
0312-XL	XL	16-17 in / 40.6-43.2 cm
0312-2XL	2XL	17-18 in / 43.2-45.7 cm
0312-3XL	3XL	18-19 in / 45.7-48.3 cm
0312-4XL	4XL	19-20 in / 48.3-50.8 cm
0312-5XL	5XL	20-22 in / 50.8-55.9 cm





MODEL: 0314



NEOPRENE KNEE WRAP WITH STABILIZER PAD HCPCS: A4466 • MINIMUM

- Opens flat for easy application over swollen or tender knees
- Hook compatible laminate allows for easy adjustment
- Multifunction stabilizer pad surrounds kneecap, provides compression and stabilization
- Non-perforated material retains body heat and increases circulation, helps maintain joint flexibility

INDICATIONS: Sprains and strains, Chondromalacia, osteoarthritis, patellofemoral tracking

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0314 UNIVERSAL

MODEL: 0315



NEOPRENE THIGH SUPPORT WITH OVAL PAD HCPCS: A9270 • MINIMUM

- Provides warmth and compression to injured thigh and hamstring areas
- Oval patch for added compression, helps reduce muscle strains and pulls
- Four-way stretch material provides even compression, plush lined interior wicks away perspiration

INDICATIONS: Nonspecific pain in the thigh area, hamstring strain or injury, mild muscle and tendon pain, heat retention and compression

MEASURE AROUND THE MIDDLE OF THE THIGH

0315-S	S	13-16 in / 33-40.6 cm
0315-M	M	16-21 in / 40.6-53.3 cm
0315-L	L	21-26 in / 53.3-66 cm
0315-XL	XL	26-32 in / 66-81.3 cm

MODEL: 0326



NEOPRENE PATELLA STABILIZER HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over all surface anatomy prominences of the joint
- Tubular pad provides extra uplift around the patellar tendon for stability, proper tracking and to help relieve minor pain
- Easy wraparound design for quick application

INDICATIONS: Sprains and strains, mild patella femoral tracking

MEASURE AROUND THE BEND OF THE KNEE

0326-S	S	13-14 in / 33-35.5 cm
0326-M	M	14.5-15.5 in / 36.8-39.3 cm
0326-L	L	16-17 in / 40.6-43.2 cm
0326-XL	XL	17.5-18.5 in / 44.4-46.9 cm

MODEL: 0212



KNEE WRAP WITH STABILIZING PAD HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over surface anatomy
- Hook compatible laminate allows for a wide range of adjustability over swelling
- Multi-function buttress surrounds patella and provides compression and stabilization

INDICATIONS: Mild support and compression, mild patella bursitis

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0212 UNIVERSAL

MODEL: 0211



THERAPEUTIC KNEE GUARD HCPCS: A4465 • MINIMUM

- Targets the physiological causes of minor knee pain to due to abnormal tracking or alignment of the patella
- As the knee moves, it applies pressure to the joint liner and patellar tendon to relieve minor pain
- Tubular pad provides uplift support around the patella for proper function

INDICATIONS: Osgood Schlatter Syndrome, patella tendonitis, anterior knee pain

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0211 UNIVERSAL

MODEL: 0216



NEOPRENE KNEE SUPPORT WITH PATELLAR PAD HCPCS: A4466 • MINIMUM

- Four-way stretch neoprene material retains body heat to increase circulation and maintains flexibility, provides even compression over the prominences of the knee
- Multi-function pad surrounds patella and provides targeted support and stabilization
- Opening over patella prevents binding

INDICATIONS: Sprains and strains, mild patella femoral tracking

MEASURE AROUND THE BEND OF THE KNEE

0216-S	S	13-14 in / 33-35.5 cm
0216-M	M	14-15 in / 35.5-38.1 cm
0216-L	L	15-16 in / 38.1-40.6 cm
0216-XL	XL	16-17 in / 40.6-43.2 cm





KNEE SUPPORTS

MODEL: 0057



CRISS-CROSS KNEE SUPPORT HCPCS: A4466 • MINIMUM

- Criss-cross design in back minimizes bunching, enhances wearing comfort
- Strong controlled stretch elastic provides effective compressive support over the knee
- Reinforced non-roll edges

INDICATIONS: Sprains and strains, general support

MEASURE AROUND THE BEND OF THE KNEE

0057-S	S	10-12.75 in / 25.432.4 cm
0057-M	M	13-15.75 in / 33-40 cm
0057-L	L	16-18.75 in / 40.6-47.6 cm
0057-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0062



SHEER ELASTIC KNEE SUPPORT HCPCS: A4466 • MINIMUM

- Strong yet sheer two-way stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Works in concert with the motion of the knee, encircles like a bandage but without the bulk

INDICATIONS: Mild sprains and strains, joint weakness due to previous injury or overexertion

MEASURE AROUND THE BEND OF THE KNEE

0062-S	S	10-12.75 in / 25.432.4 cm
0062-M	M	13-15.75 in / 33-40 cm
0062-L	L	16-18.75 in / 40.6-47.6 cm
0062-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 70/44



CONTOUR CUT KNEE SUPPORT HCPCS: A4466 • MINIMUM

- Provides support and confidence during activity, helps prevent reinjury
- Strong controlled stretch elastic provides effective compressive support over the knee
- Reinforced non-roll edges

INDICATIONS: Sprains and strains, general support

MEASURE AROUND THE BEND OF THE KNEE

0070/44-S	S	10-12.75 in / 25.432.4 cm
0070/44-M	M	13-15.75 in / 33-40 cm
0070/44-L	L	16-18.75 in / 40.6-47.6 cm
0070/44-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0072



KNEE BRACE WITH FLEXIBLE STAYS HCPCS: A4466 •• MEDIUM

- Spiral stays provide moderate medial/lateral stability
- Strong controlled stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Excellent for post-surgical recovery, increases wearer confidence for walking or standing

INDICATIONS: Sprains and strains, osteoarthritis, general support

MEASURE AROUND THE BEND OF THE KNEE

0072-S	S	10-12.75 in / 25.432.4 cm
0072-M	M	13-15.75 in / 33-40 cm
0072-L	L	16-18.75 in / 40.6-47.6 cm
0072-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0065



LIGHTWEIGHT KNEE SUPPORT WITH HINGED BARS HCPCS: L1812 ••• MAXIMUM

- Provides optimum support after moderate injury involving surrounding soft tissues
- Lightweight rigid metal hinged bars provide excellent side to side control and allow for free flexion and extension
- One-way stretch elastic sleeve utilizes a real rubber core as its basis of support

INDICATIONS: Mild to moderate medial/lateral knee instability, compression and support, post-operative recovery

MEASURE AROUND THE BEND OF THE KNEE

0065-S	S	10-12.75 in / 25.4-32.4 cm
0065-M	M	13-15.75 in / 33-40 cm
0065-L	L	16-18.75 in / 40.6-47.6 cm
0065-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0070



FIRM ELASTIC KNEE SUPPORT HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Anatomic design works in concert with the motion of the knee, encircles like a bandage but without the bulk

INDICATIONS: Mild sprains and strains, joint weakness due to previous injury or overexertion

MEASURE AROUND THE BEND OF THE KNEE

0070-S	S	10-12.75 in / 25.432.4 cm
0070-M	M	13-15.75 in / 33-40 cm
0070-L	L	16-18.75 in / 40.6-47.6 cm
0070-XL	XL	19-21.75 in / 48.3-55.2 cm

MODEL: 0074



KNEE BRACE WITH HOR-SHU SUPPORT PAD HCPCS: A4466 •• MEDIUM

- Hor-Shu support pad surrounds knee cap and provides compression and stabilization
- Strong controlled stretch elastic provides effective compressive support over the knee along with all-day wearing comfort
- Excellent for post-surgical recovery, increases wearer confidence for walking or standing

INDICATIONS: Patellar tendon sprains and strains, post-operative conditions, post-trauma discomfort

MEASURE AROUND THE BEND OF THE KNEE

0074-S	S	10-12.75 in / 25.432.4 cm
0074-M	M	13-15.75 in / 33-40 cm
0074-L	L	16-18.75 in / 40.6-47.6 cm
0074-XL	XL	19-21.75 in / 48.3-55.2 cm



ANKLE AND FOOT SUPPORTS



MODEL: 0307



NEOPRENE ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Four-way stretch perforated material provides even compression over surface anatomy
- Plush lined interior absorbs moisture and wicks away perspiration
- Fits easily inside shoes

INDICATIONS: Mild sprains and strains, edema, osteoarthritis, weakened ankles, post-cast support

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0307-S	S	6.75-7.75 in / 17.1-19.7 cm
0307-M	M	7.75-8.75 in / 19.7-22.2 cm
0307-L	L	8.75-9.75 in / 22.2-24.8 cm
0307-XL	XL	9.75-10.75 in / 24.8-27.3 cm

MODEL: 0313



NEOPRENE ANKLE WRAP HCPCS: A4466 • MINIMUM

- Opens flat for easy application over injured, swollen or tender ankles
- Hook compatible laminate allows for wide range of adjustment
- Criss-cross straps provide additional dorsiflexion assist

INDICATIONS: Sprains and strains, ankle instability, post cast removal

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0313	UNIVERSAL
------	-----------

MODEL: 0217



NEOPRENE ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Figure-8 design minimizes uncomfortable bunching during ankle motion
- Four-way stretch material provides even compression over surface anatomy
- Plush lined interior absorbs moisture and wicks away perspiration

INDICATIONS: Mild sprains and strains, edema, osteoarthritis, weakened ankles, post-cast support

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0217-S	S	6-8 in / 15.2-20.3 cm
0217-M	M	8-10 in / 20.3-25.4 cm
0217-L	L	10-12 in / 24.4-30.5 cm
0217-XL	XL	12-15 in / 30.5-38.1 cm

MODEL: 60/45



FIGURE 8 ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Figure-8 design minimizes uncomfortable bunching during ankle motion
- Elastic cotton-rayon blend fabric provides firm support, improves circulation to promote healing

INDICATIONS: Mild sprains and strains, edema, osteoarthritis, weakened ankles, post-cast support

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0060/45-S	S	7-8 in / 17.8-20.3 cm
0060/45-M	M	8.25-9 in / 21-22.9 cm
0060/45-L	L	9.25-10 in / 23.5-25.4 cm
0060/45-XL	XL	10.25-11 in / 26-27.9 cm
0060/45-VL	VL	11.25-12 in / 28.6-30.5 cm

MODEL: 0214



UNIVERSAL ANKLE STABILIZER HCPCS: A4466 • MEDIUM

- Provides more support than an ordinary elastic anklet but does not restrict motion altogether like some harness type supports
- Flexible support stays work with the sturdy fabric to provide maximum support while permitting normal activity

INDICATIONS: Ankle sprains, chronic ankle instability

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0214	UNIVERSAL
------	-----------

MODEL: 0215



BLACK ELASTIC ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Provides compression where it is needed most without binding or cutting around the edges
- Exerts compression over the soft tissues of the ankle, fitting comfortably around the bony prominences of the joint

INDICATIONS: Mild sprains and strains, weakened ankles, post-cast support

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0215-S/M	S/M	6.75-8.75 in / 17.1-22.2 cm
0215-L/XL	L/XL	8.75-10.75 in / 22.2-27.3 cm





ANKLE AND FOOT SUPPORTS

MODEL: 0060



FIRM ELASTIC ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the ankle along with all-day wearing comfort
- Anatomic design works in concert with the motion of the ankle, encircles like a bandage but without the bulk

INDICATIONS: Aching joints and muscles, edema, joint tenderness

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0060-S	S	6-8 in / 15.2-20.3 cm
0060-M	M	8-10 in / 20.3-25.4 cm
0060-L	L	10-12 in / 24.4-30.5 cm
0060-XL	XL	12-15 in / 30.5-38.1 cm

MODEL: 0063



ANKLE SUPPORT WITH SPIRAL STAYS HCPCS: A4466 • MEDIUM

- Helps reduce swelling and minimize pain while allowing full range of motion.
- Bilateral spiral stays move with the ankle and provide medial-lateral support and stability
- Excellent for light sports and exercise activity

INDICATIONS: Acute ankle sprains, chronic ankle instability, post cast removal

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0063-S	S	6.75-7.75 in / 17.1-19.7 cm
0063-M	M	7.75-8.75 in / 19.7-22.2 cm
0063-L	L	8.75-9.75 in / 22.2-24.8 cm
0063-XL	XL	9.75-10.75 in / 24.8-27.3 cm

MODEL: 0064



SHEER ELASTIC ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Sheer elastic fabric provides lightweight, firm support; helps improve circulation to promote healing
- Combines effective support with unbeatable wearing comfort
- Fits easily inside shoes
- Inconspicuous beige color

INDICATIONS: Mild sprains and strains, edema, osteoarthritis, weakened ankles, post-cast support

MEASURE AROUND ANKLE AT THE SMALLEST POINT

0064-S	S	6.75-7.75 in / 17.1-19.7 cm
0064-M	M	7.75-8.75 in / 19.7-22.2 cm
0064-L	L	8.75-9.75 in / 22.2-24.8 cm
0064-XL	XL	9.75-10.75 in / 24.8-27.3 cm



ELBOW STRAPS AND SUPPORTS

MODEL: 0301



TENNIS ELBOW STRAP WITH SUPPORT PAD HCPCS: A4466 • MEDIUM

- Excellent preventative/therapeutic device for wear during sports or construction activity
- Helps dampen vibration and shock often transmitted during periods of repetitive motion

INDICATIONS: Medial and/or lateral epicondylitis (tennis or golfer's elbow)

MEASURE AROUND THE BEND OF THE ELBOW

0301-S	S	7.75-9.75 in / 12.1-15.2 cm
0301-M	M	9.75-11 in / 15.2-18.4 cm
0301-L	L	11-13.5 in / 18.4-21.6 cm

MODEL: 0219



NEOPRENE ELBOW SUPPORT HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over surface anatomy
- Plush lined interior absorbs moisture and wicks away perspiration

INDICATIONS: Tendonitis, bursitis, arthritis, medial and/or lateral epicondylitis (tennis / golfer's elbow), mild elbow pain)

MEASURE AROUND THE BEND OF THE ELBOW

0219-S	S	9-10 in / 22.9-25.4 cm
0219-M	M	10-11 in / 25.4-27.9 cm
0219-L	L	11-12 in / 27.9-30.4 cm
0219-XL	XL	12-13 in / 30.4-33 cm



ELBOW STRAPS AND SUPPORTS



MODEL: 0053



PULLOVER ELASTIC ELBOW SUPPORT HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the elbow along with all-day wearing comfort
- Cotton-rayon blend fabric helps improve circulation to promote healing

INDICATIONS: Compressive support for weak or injured elbows

MEASURE AROUND THE BEND OF THE ELBOW

0053-S	S	8.5-10 in / 21.6-25.4 cm
0053-M	M	10-11.5 in / 25.4-29.2 cm
0053-L	L	11.5-13 in / 29.2-33 cm
0053-XL	XL	13-14.5 in / 33-36.8 cm

MODEL: 0302



ELBOW SUPPORT WITH ENCIRCLING STRAP HCPCS: A4466 • MEDIUM

- Four-way stretch material provides even compression over all the anatomical prominences of the elbow joint
- Opening over elbow helps support maintain correct position
- Encircling strap prevents painful overuse of underlying muscles
- Black (0302BL) or royal (0302RB)

INDICATIONS: Tendonitis, bursitis, arthritis, medial and/or lateral epicondylitis (tennis / golfer's elbow), mild elbow pain

MEASURE AROUND THE BEND OF THE ELBOW

0302-S	S	7.75-9.75 in / 12.1-15.2 cm
0302-M	M	9.75-11 in / 15.2-18.4 cm
0302-L	L	11-13.5 in / 18.4-21.6 cm
0302-XL	XL	13.5-15 in / 21.6-25.4 cm

MODEL: 0043



TENNIS ELBOW STRAP HCPCS: A4465 • MINIMUM

- Wrapped snugly below the elbow joint, strap inhibits the muscular action that can lead to pain and weakness
- Removable pressure pad provides focused compression for maximum support and relief

INDICATIONS: Medial and/or lateral epicondylitis (tennis or golfer's elbow)

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0043	UNIVERSAL
------	-----------

MODEL: 70/42



CONTOUR CUT ELBOW SUPPORT HCPCS: A4466 • MINIMUM

- Elastic, cotton-rayon blend fabric provides firm support; helps improve circulation to promote healing
- Absorbs perspiration for all day wearing comfort
- Contour cut design for better support and wearing comfort
- Inconspicuous neutral color

INDICATIONS: Compressive support for weak or injured elbows

MEASURE AROUND THE BEND OF THE ELBOW

0070/42-S	S	9-9.75 in / 22.9-24.8 cm
0070/42-M	M	10-10.75 in / 25.4-27.3 cm
0070/42-L	L	11-11.75 in / 27.9-29.8 cm
0070/42-XL	XL	12-13 in / 20.5-33 cm

WRIST SUPPORTS



MODEL: 0303



NEOPRENE WRIST-THUMB SUPPORT HCPCS: A4466 • MINIMUM

- Comfortably supports both wrist and thumb in opposition
- Plush lined interior absorbs moisture and wicks away perspiration

INDICATIONS: Tendonitis, mild sprains and strains, repetitive motion injury

MEASURE AROUND THE BREAK OF THE WRIST

0303-S	S	4.75-6 in / 12.1-15.2 cm
0303-M	M	6-8 in / 15.2-20.3 cm
0303-L	L	8-10 in / 20.3-25.4 cm

MODEL: 0128



NEOPRENE WRIST WRAP HCPCS: A4466 • MINIMUM

- Easy one-hand application
- Thumb hole maintains support in position
- Plush lined interior absorbs moisture and wicks away perspiration
- Royal blue neoprene

INDICATIONS: Tendonitis, mild sprains and strains, carpal tunnel syndrome

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0128	UNIVERSAL
------	-----------





WRIST SUPPORTS

MODEL: 0045



FIRM ELASTIC WRIST SUPPORT HCPCS: A4466 • MINIMUM

- Strong one-way stretch elastic provides effective compressive support over the wrist along with all-day wearing comfort
- Cotton-rayon blend fabric helps improve circulation to promote healing

INDICATIONS: Compressive support for weak or injured wrists

MEASURE AROUND THE BREAK OF THE WRIST

0045-S	S	5.5-6.5 in / 14-16.5 cm
0045-M	M	6.5-7.5 in / 16.5-19 cm
0045-L	L	7.5-8.5 in / 19-21.6 cm
0045-XL	XL	8.5-9.5 in / 21.6-24.1 cm

MODEL: 0305



NEOPRENE THUMB SPLINT HCPCS: L3923 ••• MAXIMUM

- Allows unrestricted movement in the fingers and preserves maximum function of the hand
- Contoured metal splint can be reshaped
- Contact closure fasteners permit easy, one-handed adjustment and closure

INDICATIONS: DeQuervain's tendonitis, MCP joint sprains and strains, basal joint arthritis, cumulative joint

MEASURE AROUND THE BREAK OF THE WRIST

0305-S	S	4.75-6 in / 12.1-15.2 cm
0305-M	M	6-8 in / 15.2-20.3 cm
0305-L	L	8-10 in / 20.3-25.4 cm

MODEL: 0218



NEOPRENE WRIST WRAP HCPCS: A4466 • MINIMUM

- Easy one-hand application
- Thumb hole maintains support in position
- Plush lined interior absorbs moisture and wicks away perspiration
- Black neoprene

INDICATIONS: Tendonitis, mild sprains and strains, carpal tunnel syndrome

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0218	UNIVERSAL
------	-----------

MODEL: 0047



WRAPAROUND WRIST SUPPORT HCPCS: A4466 • MINIMUM

- Provides maximum support over center of wrist and less tension around edges of support
- Thumb loop allows for easy one-hand application
- Elastic, cotton-rayon fabric absorbs perspiration for all day wearing comfort

INDICATIONS: Tendonitis, mild sprains and strains, carpal tunnel syndrome

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0047	UNIVERSAL
------	-----------

MODEL: 0304



NEOPRENE COCK-UP WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Holds the wrist in a neutral position to promote healing and reduce pain
- Perforated neoprene glove is breathable for comfortable, extended wear
- Two preshaped metal splints – one over the palm and the other over the back of the hand – are removable for additional contouring if necessary

INDICATIONS: Mild to moderate wrist sprains, strains, carpal tunnel syndrome, post-cast application

MEASURE AROUND THE BREAK OF THE WRIST

0304-XS	XS	3.75-4.75 in / 9.5-12.1 cm
0304-S	S	4.75-6 in / 12.1-15.2 cm
0304-M	M	6-7.25 in / 15.2-18.4 cm
0304-L	L	7.25-8.5 in / 18.4-21.6 cm
0304-XL	XL	8.5-10 in / 21.6-25.4 cm

MODEL: 0051



REVERSIBLE CLOTH WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Made from durable white cotton fabric – provides longer wear and more support than an ordinary elastic wristlet
- Preshaped aluminum splint keeps wrist supported at the proper angle to promote healing, relieve pain
- Back lace arrangement with contact closure fastening

INDICATIONS: Mild to moderate wrist sprains, strains, carpal tunnel syndrome, post-cast application

MEASURE AROUND THE BREAK OF THE WRIST

0051-S	S	5.5-7.5 in / 14-19 cm
0051-M	M	6.5-8.5 in / 16.5-21.6 cm
0051-L	L	7.5-9.5 in / 19-24.1 cm
0051-XL	XL	8.5-10.5 in / 21.6-26.7 cm

Specify for RIGHT or LEFT wrist



WRIST SUPPORTS



MODEL: 0050



REVERSIBLE ELASTIC WRIST SPLINT HCPCS: L3908 ●●● MAXIMUM

- Elastic, cotton-rayon tri-tension fabric places greatest amount of support directly over the wrist
- Preshaped aluminum splint keeps wrist supported at the proper angle to promote healing, relieve pain
- Contact closure encircling strap permits easy self-adjustment

INDICATIONS: Mild to moderate wrist sprains, strains, carpal tunnel syndrome, post-cast application

MEASURE AROUND THE BREAK OF THE WRIST

0050-S	S	5.5-7.5 in / 14-19 cm
0050-M	M	6.5-8.5 in / 16.5-21.6 cm
0050-L	L	7.5-9.5 in / 19-24.1 cm
0050-XL	XL	8.5-10.5 in / 21.6-26.7 cm

MODEL: 50/33



ELASTIC WRIST SPLINT HCPCS: L3908 ●●● MAXIMUM

- Specially shaped splint keeps the wrist fixed at an angle to promote healing, minimize pain
- Tri-tension fabric places greatest amount of support directly over the wrist
- Contact closures allow easy self-adjustment for excellent support and comfortable fit
- Padded closure between index finger and thumb

INDICATIONS: Mild to moderate wrist sprains, strains, carpal tunnel syndrome, post-cast application

MEASURE AROUND THE BREAK OF THE WRIST

50/33-XS	XS	3.5-5.5 in / 14-8.8 cm
50/33-S	S	5.5-6.5 in / 14-16.5 cm
50/33-M	M	6.75-7.5 in / 17.1-19 cm
50/33-L	L	7.75-8.5 in / 19.6-21.6 cm
50/33-XL	XL	8.75-10.5 in / 22.2-26.6 cm

Specify for RIGHT or LEFT wrist

MODEL: 0017



CRADLE ARM SLING

HCPCS: A4565 ●● MEDIUM

- Provides well balanced support to the arm during recuperation from fracture, sprain or surgery, evenly distributing its weight
- Reversible for right or left arm, with wide range of adjustment to easily fits most adults
- Thumb loop supports the hand, maintains the wrist in a neutral position for wearing comfort

INDICATIONS: Cast support, injury to the arm, wrist or hand

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0017	UNIVERSAL
------	-----------

SHOULDER SUPPORTS



MODEL: 0327



NEOPRENE SHOULDER SUPPORT HCPCS: A4466 ●● MEDIUM

- Protects and supports the shoulder following injury
- Retains body heat and increases circulation, helps maintain flexibility
- Includes pocket for optional use of hot/cold gel pack

INDICATIONS: Glenohumeral (shoulder) sprain and strain, tendonitis of biceps and rotator cuff, rotator cuff strain, bursitis, subluxation

MEASURE AROUND CHEST AT THE LEVEL OF THE ARM PIT

0327-S	S	UP TO 40 in / 102 cm
0327-M	M	41-51 in / 104-127 cm
0327-L	L	51-60 in / 130-152 cm

MODEL: 6147



SHOULDER BRACE HCPCS: L3650 ●● MEDIUM


- Constant reminder to sit, walk and stand with proper posture
- Padded shoulder straps pass under the arms snugly and criss-cross in back
- Allows for normal movement without binding
- Degree of support is adjustable by tightening the shoulder straps

INDICATIONS: Slumping shoulders, poor posture

MEASURE AROUND CHEST AT THE LEVEL OF THE ARM PIT

6147-S	S	UP TO 40 in / 102 cm
6147-M	M	41-51 in / 104-127 cm
6147-L	L	51-60 in / 130-152 cm





BACK AND RIB SUPPORTS

MODEL: 0007



SACRO BRACE WITH THERMOPAD

HPCPS: A4466 •• MEDIUM

- Unique construction provides optimum support and comfort to the sacrum and the lower spine
- Adjustable tapered design conforms to effectively support virtually every adult figure type, male or female
- Removable ThermoPad provides added warmth for aching muscles, joints, tendons and discs that so often accompany low back pain
- 6" deep in front, 8" deep in back

INDICATIONS: Support and compression for strains and sprains of the lower back and abdomen

UNIVERSALLY ADJUSTABLE TO FIT MOST ADULTS

0007 UNIVERSAL 32-44 in / 81.3-111.8 cm

MODEL: 6170



RIB BELT FOR WOMEN

HPCPS: A4466 • MINIMUM

- White canvas material with soft porous foam padding absorbs perspiration for comfortable all-day wear
- Sturdy elastic back panel allows for easy movement and minimizes binding
- Convenient contact closure fastening for easy self-application and adjustment

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE CHEST CIRCUMFERENCE
AT THE BASE OF THE BREAST BONE

6170-S S 23-26 in / 58.4-66 cm

6170-M M 28-32 in / 71.1-81.3 cm

6170-L L 34-38 in / 86.4-96.5 cm

6170-XL XL 40-44 in / 101.6-111.8 cm

MODEL: 6171



RIB BELT FOR MEN

HPCPS: A4466 • MINIMUM

- White canvas material with soft porous foam padding absorbs perspiration for comfortable all-day wear
- Sturdy elastic back panel allows for easy movement and minimizes binding
- Convenient contact closure fastening for easy self-application and adjustment

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

MEASURE CHEST CIRCUMFERENCE
AT THE BASE OF THE BREAST BONE

6171-S S 26-32 in / 66-81.3 cm

6171-M M 34-40 in / 86.4-101.6 cm

6171-L L 42-48 in / 106.7-121.9 cm

6171-XL XL 50-60 in / 127-142.2 cm

MODEL: 0205/0206



INDUSTRIAL BELT

HPCPS: L1499 •• MEDIUM

- Allows freedom of movement, promotes proper lifting technique
- Worn over clothing, easy to apply and remove
- Elastic and foam construction is lightweight and durable
- Comfortable and effective
- Available with shoulder straps (0205S, 0206S) or without (0205, 0206), in black (0205, 0205S) or white (0206, 0206S)

INDICATIONS: Prevention of strain and injury during heavy lifting and work activity

MEASURE WAIST CIRCUMFERENCE

		WAIST	STRAPS	COLOR
0205-XS	XS	26 - 36"	NO	BLACK
0205-R	R	34 - 44"	NO	BLACK
0205-XL	XL	42 - 52"	NO	BLACK
0205-2XL	2XL	50 - 60"	NO	BLACK
0205S-XS	XS	26 - 36"	YES	BLACK
0205S-R	R	34 - 44"	YES	BLACK
0205S-XL	XL	42 - 52"	YES	BLACK
0205S-2XL	2XL	50 - 60"	YES	BLACK
0206-XS	XS	26 - 36"	NO	WHITE
0206-R	R	34 - 44"	NO	WHITE
0206-XL	XL	42 - 52"	NO	WHITE
0206-2XL	2XL	50 - 60"	NO	WHITE
0206S-XS	XS	26 - 36"	YES	WHITE
0206S-R	R	34 - 44"	YES	WHITE
0206S-XL	XL	42 - 52"	YES	WHITE
0206S-2XL	2XL	50 - 60"	YES	WHITE



HERNIA AND SPORTS SUPPORTS



MODEL: 0005



HERNIA BELT HCPCS: L8310

•• MEDIUM

- Can be worn with both adjustable pads for a bilateral hernia or with either pad removed for a single hernia
- Complete range of pad adjustment
- Soft knit leg straps maintain pads in proper position
- Adjustable belt maintains position on the body for maximum support, containment and control

INDICATIONS: Support and control of a reducible inguinal or scrotal hernia

MEASURE HIP CIRCUMFERENCE

0005-S	S	26-32 in / 66-81.3 cm
0005-M	M	34-40 in / 86.4-101.6 cm
0005-L	L	42-48 in / 106.7-121.9 cm

MODEL: 0081



SPORTS SUPPORT

HCPCS: A9270 • MINIMUM

- Wide waist band prevents shifting on body during strenuous activity
- Soft textured pouch prevents strain and fatigue without binding
- Comfortable non-roll leg straps
- Quick dry material for easy care

INDICATIONS: Scrotal support during athletics and other strenuous activities, post-operative support

MEASURE WAIST CIRCUMFERENCE

0081-S	S	20-26 in / 50.8-66 cm
0081-M	M	26-32 in / 66-81.3 cm
0081-L	L	32-38 in / 81.3-96.5 cm
0081-XL	XL	38-44 in / 96.5-111.7 cm

MODEL: 0052



SUSPENSORY •• MEDIUM

- Soft knitted cotton/rayon mesh pouch, breathable for cool comfort and support
- Pouch holds position when body moves, is inconspicuous under clothing
- Non-elastic adjustable waist band with easy snap fastener – snaps on and off without readjustment
- Sanitary and completely washable

INDICATIONS: Relief from strain or fatigue over the scrotal area, following hernia surgery, management of scrotal disease

SIZE IS DETERMINED BY WIDTH OF POUCH AND WAIST CIRCUMFERENCE

0052-S	S	3.5-4 in / 8.9-10.1 cm UP TO 32 in / 81.3 cm WAIST
0052-M	M	4-4.5 in / 10.1-11.4 cm UP TO 36 in / 91.4 cm WAIST
0052-L	L	4.5-5 in / 11.4-12.7 cm UP TO 40 in / 101.6 cm WAIST
0052-XL	XL	5-6 in / 12.7-15.2 cm UP TO 46 in / 116.8 cm WAIST
0052-2XL	2XL	6-7 in / 15.2-17.8 cm UP TO 52 in / 132.1 cm WAIST

MODEL: 0086



SPORTS BRIEF

HCPCS: A9270 • MINIMUM

- All the protection of an athletic support with all the comfort of regular underwear
- Cool, quick drying fabrics
- Reinforced waist band
- Comfortable, non-binding leg openings
- Appropriate for all sports

INDICATIONS: Scrotal support during athletics and other strenuous activities, post-operative support

MEASURE WAIST CIRCUMFERENCE

0086-R	R	28-36 in / 71.1-91.4 cm
0086-L	L	38-44 in / 96.5-111.8 cm

HOT AND COLD COMPRESSES




MODEL: 5032



MODEL: 5034



MODEL: 5035



THERMA-KOOL COMPRESS KITS

- Reusable, microwaveable, non-toxic, latex free.
- Compress stays soft after repeated usage
- Versatile – can be used as either a hot or cold pack
- Soft and pliable even at 0° Fahrenheit
- Neat and clean – no dripping or mess
- Lightweight, reusable and economical
- Each compress comes with its own removable cover with fastening strap

INDICATIONS: Immediate post-injury and recovery phases

SELECT BASED ON APPLICATION

5032	4 x 9 in	ANKLE, ELBOW, KNEE
5034	6 x 10 in	KNEE, SHOULDER
5035	8 x 10 in	LOW BACK, SHOULDER





BANDAGES AND CAST PROTECTORS

MODEL: 0122



MODEL: 0123



MODEL: 0124



ELASTIC BANDAGES

HPCPS: A4590 • MINIMUM

- Allows for adjustable compression over injured soft tissues
- Professional quality, highly stretchable beige elastic with two fasteners
- Inconspicuous neutral color
- Can be washed and reused

INDICATIONS: Post-injury applications

SELECT BASED ON APPLICATION

0122	2" WIDE	ANKLE, FOOT, WRIST
0123	3" WIDE	ANKLE, FOOT, LEG
0124	4" WIDE	LEG, SHOULDER

MODEL: 0132



MODEL: 0133



MODEL: 0134



SELF-ADHERING ELASTIC BANDAGES

HPCPS: A4590 • MINIMUM

- Self-adhering – requires no fasteners
- Allows for adjustable compression over injured soft tissues
- Inconspicuous neutral color

INDICATIONS: Post-injury applications

SELECT BASED ON APPLICATION

0132	2" WIDE	ANKLE, FOOT, WRIST
0133	3" WIDE	ANKLE, FOOT, LEG
0134	4" WIDE	LEG, SHOULDER

MODEL: 0156



CAST PROTECTOR, HALF LEG SIZE

HPCPS: A9270

- For casts and bandages below the knee
- Water tight protection keeps cast or bandage dry while in shower or bath, or even while swimming
- Made from puncture and tear resistant vinyl
- Adjustable gripper strip provides maximum seal

MEASURE THE LENGTH OF LEG TO BE COVERED

0156-Y	Y (YOUTH)	UP TO 18 in / 45.7 cm
0156-A	A (ADULT)	UP TO 27 in / 68.6 cm

MODEL: 0157



CAST PROTECTOR, FULL LEG SIZE

HPCPS: A9270

- For casts and bandages above the knee
- Water tight protection keeps cast or bandage dry while in shower or bath, or even while swimming
- Made from puncture and tear resistant vinyl
- Adjustable gripper strip provides maximum seal

MEASURE THE LENGTH OF LEG TO BE COVERED

0157-Y	Y (YOUTH)	UP TO 24 in / 60.1 cm
0157-A	A (ADULT)	UP TO 40 in / 101.6 cm

MODEL: 0158



CAST PROTECTOR, HALF ARM SIZE

HPCPS: A9270

- For casts and bandages below the elbow
- Water tight protection keeps cast or bandage dry while in shower or bath, or even while swimming
- Made from puncture and tear resistant vinyl
- Adjustable gripper strip provides maximum seal

MEASURE THE LENGTH OF ARM TO BE COVERED

0158-Y	Y (YOUTH)	UP TO 18 in / 45.7 cm
0158-A	A (ADULT)	UP TO 24 in / 60.1 cm

MODEL: 0159



CAST PROTECTOR, FULL ARM SIZE

HPCPS: A9270

- For casts and bandages above the elbow
- Water tight protection keeps cast or bandage dry while in shower or bath, or even while swimming
- Made from puncture and tear resistant vinyl
- Adjustable gripper strip provides maximum seal

MEASURE THE LENGTH OF ARM TO BE COVERED

0159-Y	Y (YOUTH)	UP TO 24 in / 60.1 cm
0159-A	A (ADULT)	UP TO 32 in / 81.3 cm





CHAMPION KidsLine is a family of Pediatric and Youth sized braces and supports. They follow the same anatomically correct profile as adult products but are “right sized” to fit a smaller frame. Pediatric sizes will fit children from 2 to 7 years of age. Youth sizes will fit kids from 8 to 12 years old. Ages are approximate and should be used only as a quick reference guide for sizing.

MODEL: 0316



KIDS LINE

KNEE SLEEVE WITH OPEN PATELLA HCPCS: A4466 • MINIMUM

- Four-way stretch material provides even compression over the surface anatomy of the joint
- Plush lined interior absorbs moisture, wicks away perspiration

INDICATIONS: Sprains and strains, osteoarthritis

MEASURE AROUND THE BEND OF THE KNEE

0316-S	S	6-8 in
0316-M	M	8-10 in
0316-L	L	10-12 in

MODEL: 0317



KIDS LINE

ANKLE SUPPORT HCPCS: A4466 • MINIMUM

- Four-way stretch perforated material provides even compression over surface anatomy
- Fits easily inside shoes

INDICATIONS: Mild sprains and strains, post-cast support

SMALLEST PART OF THE ANKLE

0317-P	P (PEDIATRIC)	6.5-7.25 in
0317-Y	Y (YOUTH)	7.25-8.5 in

MODEL: 0318



KIDS LINE

REVERSIBLE WRIST WRAP HCPCS: A4466 • MINIMUM

- Easy one-hand application
- Thumb hole maintains support in position
- Plush lined interior absorbs moisture and wicks away perspiration

INDICATIONS: Tendonitis, mild sprains and strains, carpal tunnel syndrome

UNIVERSALLY ADJUSTABLE TO FIT MOST CHILDREN

0318	UNIVERSAL
------	-----------

MODEL: 0320



KIDS LINE

CRADLE STYLE ARM SLING HCPCS: A4565 • MINIMUM


- Provides well balanced support to the arm, evenly distributing its weight
- Thumb loop feature supports the hand, maintains the wrist in a neutral position
- Reversible for right or left arm
- Available in navy (0320NV) or print (0320PP)

INDICATIONS: For wear post-fracture, sprain or surgery

MEASURE FROM ELBOW TO TIP OF LITTLE FINGER

0320-I	I (INFANT)	3.5-5.5 in
0320-P	P (PEDIATRIC)	5-8 in
0320-Y	Y (YOUTH)	7-10 in

MODEL: 0321



KIDS LINE

FOAM CERVICAL COLLAR HCPCS: L0210 • MINIMUM

- Soft density foam – 1.5” thick – covered in stockinette
- Especially convenient for wear while sleeping
- Contact closure fasteners permit quick and easy adjustment

INDICATIONS: Treatment of minor neck pain or injuries, and as a reminder against making sudden and painful neck movements

MEASURE NECK CIRCUMFERENCE

0321-I	I (INFANT)	9-10 in
0321-P	P (PEDIATRIC)	10-12 in
0321-Y	Y (YOUTH)	12-14 in

MODEL: 0322



KIDS LINE

SOFT FOAM WRIST SPLINT HCPCS: L3908 ••• MAXIMUM

- Preshaped rigid ultra lightweight splint
- Glove with soft lined interior is breathable for comfort
- Easy to apply and remove with one hand

INDICATIONS: Mild to moderate wrist sprains and strains, post cast removal

MEASURE AROUND THE BREAK OF THE WRIST

0322-I	I (INFANT)	UNDER 4.5 in
0322-P	P (PEDIATRIC)	4.5-5.5 in
0322-Y	Y (YOUTH)	5.5-6.5 in

MODEL: 0323



KIDS LINE

ROCKER BOTTOM CAST SHOE HCPCS: L3908 ••• MAXIMUM

- Made from colorful clown print canvas with three closures
- Skid resistant rocker bottom
- Open heel and toe for versatility in sizing

INDICATIONS: To shield a lower extremity cast from moisture, shock and everyday wear


SOLE DIMENSIONS

0323-I	I (INFANT)	5.5 in IN LENGTH 3 in IN WIDTH
0323-P	P (PEDIATRIC)	7.25 in IN LENGTH 3.625 in IN WIDTH





BRACING AND SUPPORTS PRODUCT INDEX



SAI's planogram layouts make it easier for you to merchandise, for your fitters to sell from, and for your more self-service clients to shop. Modular planogram components enable you to grow your business gradually, purchasing at your own pace, and getting the most return for your investment. Choose between starter and intermediate assortments as well as brand specific assortments, depending on your business' particular needs and focus.



1791	11	2425	7	2688	6
1792	11	2426	10	2689	17
1793	10	2427	12	2786	15
2072	13	2452	18	2884	15
2081	13	2453	19	2885	6
2082	13	2454	5	2886	16
2083	13	2455	18	2890	15
2086	13	2458	6	2891	15
2087	13	2459	17	2892	15
2088	20	2460	12	2893	15
2089	12	2464	19	2894	16
2092	9	2465	19	2895	16
2095	10	2466	19	2936	18
2096	10	2501	19	2939	18
2351	14	2503	5	2956	17
2360	14	2508	6	2955	6
2362	14	2509	17	2958	18
2364	14	2518	17	2958D	18
2371	9	2540	7	5032	31
2372	9	2541	4	50/33	29
2375	9	2542	4	5034	31
2376	9	2543	4	5035	31
2382	5	2544	4	60/45	25
2383	5	2545	8	6147	29
2386	5	2546	4	6170	30
2387	5	2547	10	6171	30
2389	14	2548	7	70/42	27
2394	19	2549	7	70/44	24
2416	7	2553	8	7901Q	20
2417	9	2554	8	79030	20
2418	14	2555	8	79600	20
2419	12	2555V	8	79020	20
2421	12	2560	10	79040	20
2422	7	2658	6	CMS-001-17	11
2423	12	2659	17	CMS-001-11	11
				CMS-003-17	11
				CMS-003-11	11



0005	31	0159	32
0007	30	0205	30
0017	29	0206	30
0043	27	0211	23
0045	28	0212	23
0047	28	0214	25
0050	29	0215	25
0051	28	0216	23
0052	31	0217	25
0053	27	0218	28
0057	24	0219	26
0060	26	0301	26
0062	24	0302	27
0063	26	0303	27
0064	26	0304	28
0065	24	0305	28
0070	24	0306	21
0072	24	0307	25
0074	24	0308	22
0081	31	0309	21
0086	31	0310	22
0122	32	0311	22
0123	32	0312	22
0124	32	0313	25
0128	27	0314	23
0132	32	0315	23
0133	32	0316	33
0134	32	0317	33
0141	21	0318	33
0142	21	0320	33
0143	21	0321	33
0144	21	0322	33
0156	32	0323	33
0157	32	0326	23
0158	32	0327	29



SAI THERAPEUTIC BRANDS